IDEXX VetStat* Calibration Gas

06-05115-01

Contents: One (1) gas bottle
FOR IN-VITRO DIAGNOSTIC USE.

Intended Use: The calibration gas is intended for use only in the IDEXX VetStat Analyzer.

Product Description: One (1) 2-L gas bottle at 9.7 bars (<145 psi) at 21°C.

Active Ingredients: 6% CO₂, 14% O₂, balance N₂

Storage and Handling: The gas bottle is stable when stored at 4°–30°C (39°–86°F) until the expiration date printed on the label.

Installation:

1. Turn on the VetStat analyzer and wait for the Home screen or the gas error message to appear.
2. Remove and discard the empty gas bottle (treat it as any low-pressure consumable).
3. Record the installation date on the new gas bottle's base for future reference.
4. Remove the cap from the new gas bottle.
5. Tap **OK** at the "Stop! Low Gas. Replace Now." error message.
6. When prompted, swipe the gas bottle bar code (located on the front of this insert). The "Please Insert Gas Bottle" message appears.
7. Insert the gas bottle in its housing and turn it clockwise until it is snug. Then tap **OK**.

Warning: Do not overtighten the gas bottle as this may cause damage to the VetStat analyzer.

8. When the "New Gas Bottle?" message appears, tap **Yes**. The VetStat analyzer begins to warm up and performs a gas purge, which is indicated on the bottom of the screen. When the warm-up is complete, the Home screen appears.
9. Verify that the percentage of gas in the new gas bottle is 99% (upper left of the Home screen).

Troubleshooting:

- If the VetStat analyzer does not advance past the "Insert Gas Bottle" message, check that the gas canister is fully threaded into the housing.
- If the VetStat reports that the gas bottle is expired, the gas bottle is either at the end of its shelf life or "In-Use" life, and you need to install a new gas bottle.
- If the VetStat will not read the bar code via the bar-code reader, you can enter the bar-code number via the keypad.

Caution: Pressurized container. Protect from sunlight and do not expose to temperatures exceeding 50°C (122°F). Do not pierce or burn, even after use. Do not spray on a naked flame or any incandescent material.

Disposal: Please dispose of according to the safety regulations applicable at your facility.

Performance Characteristics: For information regarding performance characteristics, see to the VetStat Operator's Guide.

This product is warranted to perform as described in the labeling and in IDEXX literature, and IDEXX disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall IDEXX be liable for any consequential damages arising out of the aforesaid express warranty.

For technical service, contact IDEXX.

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*VetStat is a trademark of IDEXX Laboratories, Inc. in the United States and/or other countries.
© 2005 IDEXX Laboratories, Inc. All rights reserved.
*VetStat è un marchio di proprietà d'IDEXX Laboratories, Inc. e protetto negli Stati Uniti e/o in altri paesi.
© 2005 IDEXX Laboratories, Inc. Tutti i diritti sono riservati.