SNAP™ 4Dx™ Plus Test

In vitro diagnostic for the detection of *Dirofilaria immitis* antigen, antibody to *Anaplasma phagocytophilum*, antibody to *Anaplasma platis*, antibody to *Borrelia burgdorferi*, antibody to *Ehrlichia canis*, and antibody to *Ehrlichia ewingii* in canine serum, plasma, or whole blood.

- Precautions and warnings**
- All wastes should be properly decontaminated prior to disposal.
 - Do not mix components from kits with different serial numbers.
 - Do not use a SNAP device that has activated prior to the addition of sample.
 - Male-only heartworm infectious potency antigen levels that are below the detection capability of this antigen test.
 - Refer to country specific Safety Data Sheet for regional hazard identification.

- Storage**
- Store at 2-8°C until expiration date.
 - SNAP devices and reagents can be stored at room temperature (18-25°C) for 6 months or until the expiration date, whichever occurs first.
 - After SNAP devices and reagents are removed from 2-8°C for **more than 24 hours**, the expiration date is 6 months or the printed expiration date, whichever occurs first. If the 6 month expiration date occurs prior to the printed expiration date, record the new date in the space provided on the kit.

Item	Reagents	Quantity
1	1 or 5 bottles Anti-D. <i>immitis</i> /A. <i>phagocytophilum</i> /A. <i>platis</i> /B. <i>burgdorferi</i> /E. <i>canis</i> /E. <i>ewingii</i> , HRPPO conjugate (preserved with ProClin™ 150)	7.0 mL
2	SNAP device	5, 15, 30, or 180

Reagents contained in each device:

Wash solution (preserved with ProClin™ 150) 0.4 mL
Substrate solution 0.6 mL


Other components: transfer pipettes, sample tubes, and reagent rack

Conjugate - H317/H412/P261/P272/P273/P280/P302+P330/P333+P313/P362+P364/P501: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapors. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves if on skin. Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. Dispose of contents/container in accordance with local/regional/national/international regulations.


Wash Solution - H317/H412/P261/P272/P273/P280/P302+P330/P333+P313/P362+P364/P501: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapors. Wash thoroughly after handling. Contaminated work clothing must not be allowed out of the workplace. Avoid release to the environment. Wear eye protection/face protection. Wear protective gloves if on skin. Wash with plenty of water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents/container in accordance with local/regional/national/international regulations.

- Sample information.**
- Samples must be at room temperature (18-25°C) before beginning the test procedure.
 - Serum, plasma, or anti-coagulated whole blood (e.g., EDTA, heparin), either fresh or stored at 2-8°C for up to one week, can be used.
 - For longer storage, serum, or plasma can be frozen (-20°C or colder) and then re centrifuged before use.
 - Hemolyzed or lipemic samples will not affect test results.

- Test procedure**
- If stored in a refrigerator, allow all components to equilibrate at room temperature (18-25°C) for 30 minutes before use. **Do not heat.**

The sample will flow across the result window, reaching the activation circle in 30-60 seconds. Some sample may remain in the sample well.


Note: Some samples may not flow to the activation circle within 60 seconds, and therefore, the circle may not turn color. In this case, press the activator after the sample has flowed across the result window.

- Read test result at **6 minutes**.
- Note:** The positive control may develop sooner, but results are not complete until 6 minutes.

Interpreting the test results

Positive result


Any color development in the sample spots indicates the presence of *Dirofilaria immitis* antigen, *A. phagocytophilum* antibody, *A. platis* antibody, *B. burgdorferi* antibody, *E. canis* antibody, or *E. ewingii* antibody in the sample.

- Notes:**
- The *A. phagocytophilum*/A. *platis* spot cannot differentiate between the two species; a positive result indicates presence of antibodies to *A. phagocytophilum* and/or *A. platis*.
 - The *E. canis*/E. *ewingii* spot cannot differentiate between the two species; a positive result indicates presence of antibodies to *E. canis* and/or *E. ewingii*.
 - In a low percentage of samples (0.027% as reported), interfering substances in the patient's blood can cause all spots on the device to react as positive. In this event, the sample should be re-assayed as plasma or serum to reduce the likelihood of interference.

Negative result

Only positive control spot develops color.

Invalid results



- Background**—If the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No color development**—If the positive control does not develop color, repeat the test.

Borrelia burgdorferi vaccine cross-reactivity—The *Borrelia burgdorferi* assay detects antibodies induced as a result of natural infection with the organism and not following immunization with the following vaccines: Recombitek™ Lyme, LymeVax™, and Galaxy™ Lyme and Novivac™ Lyme. Recombitek is a registered trademark of Merial, Inc. LymeVax is a registered trademark of Pfizer, Inc. Galaxy and Novivac are registered trademarks of Merck Animal Health.

Symbol	Descriptions
Use by date	Authorized representative in the European Community
Batch code (lot)	Consult instructions for use
Serial number	In vitro diagnostic
Temperature limitation	Diagnostic in vitro
Manufacturer	Date of manufacture
Catalog number	

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