


Diagnostic algorithm for diagnosis and staging of FeLV infection—Europe



For cats used as a blood donors or for breeding, screening with FeLV Quant RealPCR Test and FeLV Antigen by ELISA is recommended.

*Lymphoma and bone marrow suppression has been described occasionally in recessively infected cats. The FeLV Quant RealPCR Test can be considered.


*Early infection can lead to a negative result for both the in-clinic and reference laboratory FeLV antigen testing options; recheck in 6 weeks if cat recently acquired or at risk of recent exposure.

^{*}Feline Triple, RealPCR and SNAP are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.



Follow-up testing (after 6 weeks) with FeLV Antigen by ELISA and FeLV Quant RealPCR* Test—Europe

Staging of infection based on test results



*Cats can shift between progressive and regressive infection over the course of the disease. If clinical signs develop, retesting with quantitative real-time PCR is recommended to confirm evidence of progression versus comorbidity with another disease.

§Proviral DNA load (copies/mL) according to Beall et al.12

References

- 1. Beall MJ, Buch J, Cahill RJ, et al. Evaluation of a quantitative enzyme-linked immunosorbent assay for feline leukemia virus p27 antigen and comparison to proviral DNA loads by real-time polymerase chain reaction. Comp Immunol Microbiol Infect Dis. 2019;67:101348. doi:10.1016/j.cimid.2019.101348
- 2. Beall MJ, Buch J, Clark G, et al. Feline leukemia virus p27 antigen concentration and proviral DNA load are associated with survival in naturally infected cats. Viruses

© 2023 IDEXX Laboratories, Inc. All rights reserved. • 09-2690546-00 • The IDEXX Privacy Policy is available at idexx.com.

*RealPCR is a trademark or registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

PCR testing is a service performed pursuant to an agreement with Roche Molecular Systems, Inc.

 $^{^{\}dagger}$ Abortive cats are FeLV antibody positive; IDEXX Reference Laboratories does not offer an FeLV antibody test.