




AKI to CKD and everything in between: getting to grips with confusing renal diagnostics

Nancy A. Sanders, DVM
DACVIM(SAIM), DACVECC
IDEXX Medical Education Specialist



Financial Disclosure

I am a full-time IDEXX; The relationship does not influence the medicine. It *may* influence a few of my images...


The information contained herein is intended to provide general guidance only. Diagnosis, treatment, and monitoring should be patient specific and is the responsibility of the veterinarian providing primary care.




Learning objectives

By the end of this presentation, participants should be able to:

1. Compare and contrast the significance of kidney injury vs. chronic kidney disease
2. Identify biomarkers of acute and active kidney injury, including urinary cystatin B
3. Understand the cellular sources for cystatin B and the significance of elevated urinary concentrations
4. Apply acute injury concepts to clinical cases



Timeline of renal biomarkers



Acute kidney injury (AKI) vs. chronic kidney disease (CKD): Why do we care?

AKI

- + Early detection to prevent progression
- + Institute supportive care and specific therapy when possible
- + Determine resolution or progression
- + Short-term financial and emotional investments are intense
- + Prolonged hospitalization: associated with higher morbidity and mortality

CKD

- + Early detection/intervention in attempt to slow progression
- + Institute dietary therapy, supportive care
- + Determine likelihood of rapid progression
- + Long-term financial, emotional, and time commitments
- + Usually outpatient therapy; when hospitalization is required, usually associated with low morbidity and mortality

Hallmarks of AKI (vs. CKD)

History and physical exam

- + Acute onset—hours to days
- + Toxin exposure (lily, grapes, NSAIDs, anesthetics...)
- + Renomegaly, renal pain
- + Lack of other PE change
- + Bradycardia (if severe hyperkalemia)
- + Hypothermia


Lab findings

- + Hyperkalemia
- + Urinary granular casts, normoglycemic glucosuria

Imaging

- + Renomegaly in 70%
- + Hydroureter, pyelectasia, hydronephrosis
- + Ureteral calculi
- + Normal parathyroid gland

Performance of current renal “functional” biomarkers (estimates of GFR)




Source: Hall JA, Yerramilli M, Obare E, Yerramilli M, Almes K, Jewell DE. Serum concentrations of symmetric dimethylarginine and creatinine in dogs with naturally occurring chronic kidney disease. *J Vet Intern Med.* 2016;30(3):794–802. doi:10.1111/jvim.13942

GFR biomarkers fall short
as early detectors of kidney
disease



Categorization of biomarkers and analytes used to evaluate kidney function and injury

| Indirect markers of function | Urine-based markers | Other important analytes | Acute kidney injury markers |
|---|---|---|--|
| <p>Most specific (limited extrarenal influences):</p> <ul style="list-style-type: none"> + SDMA + Creatinine <p>Less specific (more extrarenal influences):</p> <ul style="list-style-type: none"> + BUN + Phosphorus | <p>Urinalysis</p> <ul style="list-style-type: none"> + Physical + Chemical + Sediment <p>UPC</p> | <ul style="list-style-type: none"> + Potassium + Sodium/chloride + Calcium + Albumin/TP + Hematocrit + FGF-23 | <div style="border: 2px solid purple; padding: 5px;"> <ul style="list-style-type: none"> + Cystatin B + Urine Clusterin + NGAL </div> |



You need broad assessment to understand kidney health

What can we measure in clinical practice?

Glomerular function

How well are the kidneys clearing waste from the body (GFR)

Creatinine, SDMA, BUN


Tubular function

Important in solute and water management

Urine concentration and protein; serum and urine electrolytes, glucose, acid-base

Urine composition

Concentration/volume, pH, cellular and crystalline elements, infectious organisms, protein




The functional nephron

Renal tubules are where the action **really** is

Tubular function


- + The actual work of the kidney primarily takes place here. Filtering, reabsorbing, and secreting solutes and water
- + Impact urine concentration and what is excreted
- + Dysfunction can impact electrolytes, protein levels, glucose, acid-base balance
- + Captured in chemistry panel and urinalysis



Traditional renal **injury** markers are good, not great

- + Proteinuria
- + Hematuria, pyuria
- + Bacteriuria
- + Renal epithelial cells in the urine
- + Glucosuria (normoglycemia)
- + Cylindruria (casts)
- + Decreased USG

Granular casts



Cellular cast

Source: IDEXX SediVue Dx® images

Functional markers are
in blood

Tubular injury markers are
in the urine

Take-home message:

You can't assess kidney health
without urine



Veterinary criteria – IRIS AKI grading

Table 1: IRIS AKI Grading Criteria

| AKI Grade | Blood Creatinine | Clinical Description |
|------------------|--|---|
| Grade I | <1.6 mg/dl (<140 µmol/l) | <p>Nonazotemic AKI:</p> <p>a. Documented AKI: (historical, clinical, laboratory, or imaging evidence of AKI, clinical oliguria/anuria, volume responsiveness‡) and/or</p> <p>b. Progressive nonazotemic increase in blood creatinine: ≥ 0.3 mg/dl (≥ 26.4 µmol/l) within 48 h</p> <p>c. Measured oliguria (<1 ml/kg/h)# or anuria over 6 h</p> |
| Grade II | 1.7 – 2.5 mg/dl (141 – 220 µmol/l) | <p>Mild AKI:</p> <p>a. Documented AKI and static or progressive azotemia</p> <p>b. Progressive azotemic: increase in blood creatinine; ≥ 0.3 mg/dl ≥ 26.4 µmol/l) within 48 h), or volume responsiveness‡</p> <p>c. Measured oliguria (<1 ml/kg/h)# or anuria over 6 h</p> |
| Grade III | 2.6 – 5.0 mg/dl (221 – 439 µmol/l) | <p>Moderate to Severe AKI:</p> <p>a. Documented AKI and increasing severities of azotemia and functional renal failure</p> |
| Grade IV | 5.1 – 10.0 mg/dl (440 – 880 µmol/l) | |
| Grade V | >10.0 mg/dl (>880 µmol/l) | |

(‡Volume responsive is an increase in urine production to >1 ml/kg/h over 6 h; and/or decrease in serum creatinine to baseline over 48 h)

Risk
Nonazotemic

Injury
Mildly azotemic

Failure
Moderately to severely azotemic



<http://www.iris-kidney.com/education/index.html>

Subgrade

- Each grade of AKI is further subgraded as:
1. Non oliguric (NO) or oligo-anuric (O)
 2. Requiring renal replacement therapy (RRT)

<http://www.iris-kidney.com/guidelines/grading.html>

Can we do any better?


YES!!! We can!

Cystatin B bridges the gap in
our abilities to detect early
and active renal injury





What is cystatin B?


- + Member of cystatin family
 - + Protease inhibitors that help protect against leakage of proteolytic enzymes from lysosomes
 - + Trace amounts in the serum of healthy subjects
- + A small, intracellular protein
 - + 11 kDa (11,000 daltons)
 - + Ubiquitous in many cells, including proximal renal tubular cells
- + Freely filtered at the glomerulus
- + Increased urinary [cystatin B] indicates active, ongoing tubular injury
 - + Think of it as the ALT of the kidney




Cystatin B is a very small protein contained in epithelial cells of the renal tubules




During **active or acute kidney injury**, renal tubular epithelial cells can be damaged




Kidney injury markers are additive to current indirect functional markers




Kidney injury markers are additive to current indirect functional markers




Traditional diagnostics only allow for case evaluation by functional markers



Addition of an injury marker provides better case discrimination and management



Pathophysiology of AKI



Common etiologies:

- + Infectious diseases
- + Nephrotoxins
- + Systemic diseases with secondary renal involvement (inflammation)
- + Alterations in hemodynamics
- + Obstructive disorders

Causes of AKI include:

Cat

- + Toxins (plants, chemotherapeutics)
- + Pyelonephritis
- + Acute pancreatitis
- + Marked dehydration
- + Obstructive disorders
- + Etiology unknown ~30%

Dog

- + Toxins (plants, chemotherapeutics, foods)
- + Pyelonephritis
- + Acute pancreatitis
- + Marked dehydration
- + Obstructive disorders
- + Leptospirosis
- + Lyme nephritis
- + Congestive heart failure

AKI can develop in hospitalized patients: Monitor and grade daily

- + Dehydration
- + Age > very young or old
- + Diuretic or nephrotoxic drug therapy
- + Hypokalemia or hypercalcemia
- + Sepsis
- + Congestive heart failure
- + Acute pancreatitis
- + Systemic hypertension
- + CKD

Avoid iatrogenic AKI!

Nephrotoxic drugs
Hemodynamic instability
Fluid overload

**Abnormal renal functional markers reflect progression of injury to dysfunction after the fact:
Early recognition of renal injury is an opportunity to change course of disease**

Wellness labs



SDMA 14 0 – 14 µg/dL

**Stable
functional renal
markers**

6 months later

SDMA 11 0 – 14 µg/dL

**Progressive
functional
renal markers**




SDMA 14 0 – 14 µg/dL

SDMA **23** 0 – 14 µg/dL

Two dogs seen for general check up

Initial treatment of AKI



Initial therapy for AKI; back to the basics

- + Address hydration and volemic status
- + Institute disease-specific therapy whenever possible
- + Address complications of AKI
- + Avoid and monitor for complications of therapy
- + Keep the patient comfortable
- + Prevent further renal injury
- + Address nutrition

Use serial evaluations of IRIS AKI grading for objective measure of response to therapy; adapt/adjust accordingly

Fluid therapy for kidney disease: Less may be more

- + Fluids are drugs—avoid overdose
- + Fluids do not improve kidney function
- + Hypervolemia causes AKI and kills patients that already have it
- + Not every patient with kidney disease (acute or chronic) needs fluids!!!




Assessment of response to fluid therapy is essential

- + Perfusion parameters: HR, CRT, mucous membranes, pulses, lactate, base excess
- + Body weight 2-4x/day: >5%–10% increase slow or stop fluids
- + Lung auscultation: \geq q12 hrs, more frequently if any changes in RR/RE



If azotemia worsens with IV fluid therapy, consider **decreasing** fluid rate.

Especially if total daily volume exceeds maintenance or if weight gain.



My dog ate some raisins.



| Chemistry | | 3/14/24 3:32 AM | |
|---------------------------------|----------------|--------------------|--|
| Glucose | 105 | 63 - 114 mg/dL | |
| IDEXX SDMA | 10 | 0 - 14 µg/dL | |
| Creatinine | 1.0 | 0.5 - 1.5 mg/dL | |
| IDEXX SDMA | 10 | 0 - 14 µg/dL | |
| Creatinine | 1.0 | 0.5 - 1.5 mg/dL | |
| BUN | 18 | 9 - 31 mg/dL | |
| IDEXX Cystatin B (Urine) | >500 | 0 - 99 ng/mL | |
| Potassium | 5.1 | 4.0 - 5.4 mmol/L | |
| Na: K Ratio | 29 | 28 - 37 | |
| Chloride | 114 | 108 - 119 mmol/L | |
| TCO2 (Bicarbonate) | 23 | 13 - 27 mmol/L | |
| Anion Gap | 17 | 11 - 26 mmol/L | |
| Total Protein | 5.5 | 5.5 - 7.5 g/dL | |
| Albumin | 3.0 | 2.7 - 3.9 g/dL | |

3/14/24
3:32 AM

| | | |
|---------------------------|-----------------------------|-----------|
| Collection | FREECATCH | |
| Color | DARK YELLOW | |
| Clarity | TURBID | |
| Specific Gravity | 1.049 | >= 1.030 |
| pH | 5.5 | 6.0 - 7.5 |
| Urine Protein | 2+ | |
| Glucose | NEGATIVE | |
| Ketones | NEGATIVE | |
| Blood / Hemoglobin | 3+ | |
| Bilirubin | 1+ | |
| Urobilinogen | NORMAL | |
| White Blood Cells | 0-2 | |
| Red Blood Cells | 10-15 | |
| Bacteria | RARE COCCI <9/HPF | |
| Additional Bacteria | RARE RODS <9/HPF | |
| Casts | 4+ (>10)/HPF | |

Rover

- + 3-year-old MC MixB
- + Confirmed raisin ingestion
- + Amount/time prior to presentation uncertain
- + Previously healthy

Elevated urinary cystatin B, proteinuria, cylindruria

Three days later, after IV fluids for 48 hours

| Chemistry | | 3/17/24 1:07 AM | |
|--------------------------|-----|--------------------|--|
| IDEXX SDMA | 8 | 0 - 14 µg/dL | |
| Creatinine | 1.2 | 0.5 - 1.5 mg/dL | |
| BUN | 26 | 9 - 31 mg/dL | |
| IDEXX Cystatin B (Urine) | <50 | 0 - 99 ng/mL | |
| Ratio | | | |
| Phosphorus | 5.1 | 2.5 - 6.1 mg/dL | |
| Calcium | 9.5 | 8.4 - 11.8 mg/dL | |
| Sodium | 148 | 142 - 152 mmol/L | |
| Potassium | 5.1 | 4.0 - 5.4 mmol/L | |
| Na: K Ratio | 29 | 28 - 37 | |
| Chloride | 114 | 108 - 119 mmol/L | |
| TCO2 (Bicarbonate) | 25 | 13 - 27 mmol/L | |
| Anion Gap | 14 | 11 - 26 mmol/L | |
| Total Protein | 5.1 | 5.5 - 7.5 g/dL | |
| Albumin | 2.7 | 2.7 - 3.9 g/dL | |
| Globulin | 2.4 | 2.4 - 4.0 g/dL | |

| Urinalysis | | 3/17/24 1:07 AM | |
|---|--------------|--------------------|--|
| Collection | FREECATCH | | |
| J Vet Diagn Invest 17:223-231 (2005) | | | |
| Canine renal pathology associated with grape or raisin ingestion: 10 cases | | | |
| Glucose | NEGATIVE | | |
| Ketones | a RACE | | |
| Blood / Hemoglobin | 3+ | | |
| Bilirubin | 1+ | | |
| Urobilinogen | NORMAL | | |
| White Blood Cells | 0-2 | | |
| Red Blood Cells | 30-50 | | |
| Bacteria | NONE SEEN | | |
| Additional Bacteria | | | |
| Epithelial Cells | 1+ (1-2)/HPF | | |
| Mucus | NONE SEEN | | |
| Casts | NONE SEEN | | |
| Crystals | NONE SEEN | | |


Canine renal pathology associated with grape or raisin ingestion: 10 cases

All dogs had degeneration or necrosis (or both) of proximal renal tubules with basement membranes remaining intact, and epithelial regeneration was observed in 5 out of 10 cases.

Cystatin B has value with evaluating patients with CKD as well!

- + CKD progressive and irreversible
- + Rate of progression unpredictable
- + Cystatin B identifies active, progressive injury in dogs with CKD
- + Increased urinary cystatin B (uCysB) in dogs with IRIS Stage 1 CKD predicts rapid progression
- + Identifies which dogs need more frequent monitoring

Segev, et al. J Vet Intern Med. 2023




y-intercepts calculated from inverse urinary cystatin B (uCysB) vs. time

Uses for the IDEXX Cystatin B Test




Back in the day...



More contemporary view...

Your AKI patient may have or develop CKD




Your CKD patient may have concurrent active kidney injury

Kidney function in health and disease is impacted by risk factors, injury, and outcomes

Risk factors

- + Breed
- + Age
- + Sex
- + Diet
- + Drugs
- + Pre-existing disease
 - + CKD
 - + Hypertension
 - + Metabolic disease
 - + Cardiac disease



Outcomes

- + Persistent damage
- + Cardiovascular events
- + Additional kidney events
- + Diminished quality and quantity of life
- + Cost events

STAGING CKD: IRIS guidelines

| | |  |  |  |  |
|---|--------|---|---|---|---|
| | | Stage 1 No azotemia (Normal creatinine) | Stage 2 Mild azotemia (Normal or mildly elevated creatinine) | Stage 3 Moderate azotemia | Stage 4 Severe azotemia |
| Creatinine in mg/dL | | Less than 1.4 (125 µmol/L) | 1.4–2.8 (125–250 µmol/L) | 2.9–5.0 (251–440 µmol/L) | Greater than 5.0 (440 µmol/L) |
| Stage based on stable creatinine | Canine | Less than 1.4 (125 µmol/L) | 1.4–2.8 (125–250 µmol/L) | 2.9–5.0 (251–440 µmol/L) | Greater than 5.0 (440 µmol/L) |
| | Feline | Less than 1.6 (140 µmol/L) | 1.6–2.8 (140–250 µmol/L) | 2.9–5.0 (251–440 µmol/L) | Greater than 5.0 (440 µmol/L) |
| SDMA in µg/dL | | Less than 18 | 18–35 | 36–54 | Greater than 54 |
| Substage based on stable SDMA | Canine | Less than 18 | 18–35 | 36–54 | Greater than 54 |
| | Feline | Less than 18 | 18–25 | 26–38 | Greater than 38 |
| UPC ratio | | Nonproteinuric < 0.2 Borderline proteinuric 0.2–0.5 Proteinuric > 0.5 | | | |
| Substage based on proteinuria | Canine | Nonproteinuric < 0.2 Borderline proteinuric 0.2–0.5 Proteinuric > 0.5 | | | |
| | Feline | Nonproteinuric < 0.2 Borderline proteinuric 0.2–0.4 Proteinuric > 0.4 | | | |
| Systolic blood pressure in mm Hg | | Normotensive < 140 Prehypertensive 140–159 | | | |
| Substage based on blood pressure | | Hypertensive 160–179 Severely hypertensive ≥ 180 | | | |


1

2

3

Source: *Diagnosing, Staging, and Treating Chronic Kidney Disease in Dogs and Cats*. International Renal Interest Society. Accessed January 7, 2026. https://static1.squarespace.com/static/666b9ecb4064a156963b4162/t/66a6dbe460ba38633de536c0/1722211300466/IRIS_Pocket_Guide_to_CKD_2023.pdf

IRIS guidelines for staging CKD




Early marker of
chronic kidney
disease?





Yes! Fibroblast growth factor-23 (FGF-23) is an early marker of CKD

- FGF-23 is a phosphatonin; a protein/hormone that regulates phosphate balance
- Released from osteocytes in response to increased phosphorus, calcitriol, and PTH concentrations
- FGF-23 decreases phosphorus and calcitriol concentrations
- Important in animals with chronic kidney disease
- FGF-23 often rises BEFORE serum phosphorus does
- Clinical measurement is now commercially available for and indicated for cats with kidney disease (impaired phosphate excretion)



FGF-23 identifies phosphorus overload before overt hyperphosphatemia



IDEXX FGF-23 bottom line:

Persistent Findings of CKD – IRIS Stage 1 or 2 Diagnosis

Metabolic bone disease is a complex syndrome that involves fibroblast growth factor-23 (FGF-23), parathyroid hormone (PTH), 1,25-dihydroxy D3 (1,25 vitamin D₃, calcitriol), calcium, and phosphorus


FGF-23 is released in response to increases in phosphorus and calcitriol

FGF-23 often rises before total serum phosphorus

Alerting to earlier need for phosphorus reduction



Simplified medical workflow for FGF-23 in cats with CKD





Impact of renal diet on cats with early CKD: Big Data Study




Does early diagnosis and treatment of CKD with renal diet **slow progression** of disease in cats?



Does early diagnosis and treatment of CKD with renal diet **impact survival** in cats?

CKD cats treated with renal diet in early CKD stages progress more slowly




On average, early-stage CKD cats receiving renal diet remained in their current stage 11–19 months longer than similarly staged cats not fed a renal diet.

Source: Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA: IDEXX manuscript; JAVMA publication pending.

CKD cats lived longer on renal diet

Impact of renal diet on survival of CKD cats (IRIS CKD Stage 1 and 2)



Hazard ratio: **0.70** (95% CI: 0.67–0.74) $P < .001$)

Cats diagnosed and treated with renal diet in IRIS CKD Stage 1 and 2 survived 5 months longer than CKD cats not fed a renal diet.

This impact equates to an average 20% extension in survival during the 3-year study follow-up period.

Summary: Diagnosis and treatment of IRIS CKD Stage 1 and 2 cats with renal diet was associated with **slower progression** and **longer survival**.

Source: Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA: IDEXX manuscript; JAVMA publication pending.

*Increased survival time is calculated using restricted mean survival time (RMST). RMST is the area under the survival curve up to a specific time point and is generally more reliably estimable than mean or median survival times. All cause mortality was measured for survival analyses across a 3-year follow-up period.

Renal biomarkers, use every tool you have!



- **Think about limitations:**
 - Creatinine: muscle mass, early disease
 - SDMA: rare comorbidities
 - Tubular injury markers: the IDEXX Cystatin B Test
 - FGF-23: comorbidities
 - USG: comorbidities and medications
 - Proteinuria: comorbidities

Even with limitations, often by using all of the available diagnostics you can put the puzzle together.

Take Home

- Acute kidney injury and chronic kidney disease are a continuum
- A COMPLETE urinalysis is of UTMOST importance when evaluating kidney disease
- Become familiar with renal biomarkers and their indications as well as their limitations
- IDEALLY, patients at risk for renal injury or with early CKD are identified and managed BEFORE azotemia develops
- Fluid therapy paradigms have changed...dramatically.
- Newer biomarkers
 - SDMA earlier and more accurate estimate measure of GFR
 - Urine cystatin B, a urine renal biomarker, is a specific marker of ACTIVE renal tubular injury
 - Fibroblast growth factor-23 (FGF-23) is a marker of metabolic bone disease; it is an EARLY marker of CKD, and it is currently available through IDEXX for feline patients



Thank You!!



Don't just keep up.
Stay ahead with
The Curated Mind,
IDEXX's monthly
education newsletter.

THE CURATED MIND

CONTINUING
EDUCATION BROUGHT
TO YOU BY IDEXX



LEARN EVOLVE REPEAT

Scan the code to stay in the know and receive updates about:

- + Upcoming live events
- + New on-demand education resources
- + The latest *Shake Up Your Workup* podcast
- + Articles from *The Vetiverse*
- + Upcoming in-person education opportunities

