



JANUARY 17-21 • NAVC.COM • ORLANDO, FL

# WORLD CLASSIC

CELEBRATING THE CHAMPIONS OF CARE

**VMX**  
VETERINARY MEETING & EXPO



# AKI to CKD and everything in between: getting to grips with confusing renal diagnostics

Nancy Sanders DVM, DACVIM (SAIM), DACVECC


Greg Grauer DVM, MS Diplomate, ACVIM (SAIM)




# Learning objectives

**By the end of this presentation, participants should be able to:**


1. Compare and contrast the significance of kidney injury vs. chronic kidney function
2. Identify biomarkers of acute and active kidney injury, including urinary cystatin B
3. Understand the cellular sources for cystatin B and the significance of elevated urinary concentrations
4. Apply acute injury concepts to clinical cases



# Timeline of renal biomarkers



# Terminology can be confusing; ARF, AKI, ATI, ATN



J. Himmelfarb, T.A. Ikizler, Acute kidney injury: changing lexicography, definitions, and epidemiology, Kidney International, Volume 71, Issue 10, 2007, Pages 971-976, ISSN 0085-2538, <https://doi.org/10.1038/sj.ki.5002224>.  
Kellum, John A., Claudio Ronco, and Rinaldo Bellomo. "Conceptual advances and evolving terminology in acute kidney disease." Nature Reviews Nephrology 17.7 (2021): 493-502.



# Acute kidney injury (AKI) vs. chronic kidney disease (CKD): Why do we care?

## AKI

- + Early detection to prevent progression
- + Institute supportive care and specific therapy when possible
- + Determine resolution or progression
- + Short-term financial and emotional investments are intense
- + Prolonged hospitalization: associated with higher morbidity and mortality

## CKD

- + Early detection/intervention in attempt to slow progression
- + Institute dietary therapy, supportive care
- + Determine likelihood of rapid progression
- + Long-term financial, emotional, and time commitments
- + Usually outpatient therapy; when hospitalization is required, usually associated with low morbidity and mortality

# Veterinary criteria – IRIS AKI grading

**Table 1: IRIS AKI Grading Criteria**

AKI Grade	Blood Creatinine	Clinical Description
<b>Grade I</b>	<1.6 mg/dl (<140 µmol/l)	Nonazotemic AKI: a. Documented AKI: (historical, clinical, laboratory, or imaging evidence of AKI, clinical oliguria/anuria, volume responsiveness‡) and/or b. Progressive nonazotemic increase in blood creatinine: ≥ 0.3 mg/dl (≥ 26.4 µmol/l) within 48 h c. Measured oliguria (<1 ml/kg/h)# or anuria over 6 h
<b>Grade II</b>	1.7 – 2.5 mg/dl (141 – 220 µmol/l)	Mild AKI: a. Documented AKI and static or progressive azotemia b. Progressive azotemic: increase in blood creatinine; ≥ 0.3 mg/dl (≥ 26.4 µmol/l) within 48 h), or volume responsiveness‡ c. Measured oliguria (<1 ml/kg/h)# or anuria over 6 h
<b>Grade III</b>	2.6 – 5.0 mg/dl (221 – 439 µmol/l)	Moderate to Severe AKI: a. Documented AKI and increasing severities of azotemia and functional renal failure
<b>Grade IV</b>	5.1 – 10.0 mg/dl (440 – 880 µmol/l)	
<b>Grade V</b>	>10.0 mg/dl (>880 µmol/l)	

(‡Volume responsive is an increase in urine production to >1 ml/kg/h over 6 h; and/or decrease in serum creatinine to baseline over 48 h)

**Risk**  
Nonazotemic

**Injury**  
Mildly azotemic

**Failure**  
Moderately to severely azotemic



<http://www.iris-kidney.com/education/index.html>

## Subgrade

Each grade of AKI is further subgraded as:

1. Non oliguric (NO) or oligo-anuric (O)
2. Requiring renal replacement therapy (RRT)

# Hallmarks of AKI (vs. CKD)

## History and physical exam

- + Acute onset—hours to days
- + Toxin exposure (lily, grapes, NSAIDs, anesthetics...)
- + Renomegaly, renal pain
- + Lack of other PE change
- + Bradycardia (if severe hyperkalemia)
- + Hypothermia

## Lab findings

- + Hyperkalemia
- + Urinary granular casts, normoglycemic glucosuria

## Imaging

- + Renomegaly in 70%
- + Hydroureter, pyelectasia, hydronephrosis
- + Ureteral calculi
- + Normal parathyroid gland




# Back in the day...



More contemporary view...

Your AKI patient may have or develop CKD




Your CKD patient may have concurrent active kidney injury

# Kidney function in health and disease is impacted by risk factors, injury, and outcomes

## Risk factors

- + Breed
- + Age
- + Sex
- + Diet
- + Drugs
- + Pre-existing disease
  - + CKD
  - + Hypertension
  - + Metabolic disease
  - + Cardiac disease




## Potential outcomes

- + Persistent damage
- + Cardiovascular events
- + Additional kidney events
- + Diminished quality and quantity of life
- + Cost events

Source: Chawla LS, Eggers PW, Star RA, Kimmel PL. Acute kidney injury and chronic kidney disease as interconnected syndromes. N Engl J Med. 2014;371(1):58–66. doi:10.1056/NEJMra1214243

# Performance of current renal “functional” biomarkers (estimates of GFR)



GFR biomarkers fall short  
as early detectors of  
kidney disease





# Categorization of biomarkers and analytes used to evaluate kidney function and injury

Indirect markers of function	Urine-based markers	Other important analytes	Acute kidney injury markers
<p><b>Most specific (limited extrarenal impact):</b></p> <ul style="list-style-type: none"><li>+ SDMA</li><li>+ Creatinine</li></ul> <p><b>Less specific (more extrarenal impact):</b></p> <ul style="list-style-type: none"><li>+ BUN</li><li>+ Phosphorus</li></ul>	<p><b>Urinalysis</b></p> <ul style="list-style-type: none"><li>+ Physical</li><li>+ Chemical</li><li>+ Sediment</li></ul> <p><b>UPC</b></p>	<ul style="list-style-type: none"><li>+ Potassium</li><li>+ Sodium/chloride</li><li>+ Calcium</li><li>+ Albumin/TP</li><li>+ Hematocrit</li><li>+ FGF-23</li></ul>	<div><ul style="list-style-type: none"><li>+ Cystatin B</li><li>+ Urine Clusterin</li><li>+ NGAL</li></ul></div>

You need broad assessment to understand kidney health

# What can we measure in clinical practice?

## Glomerular function

How well are the kidneys clearing waste from the body (GFR)

Creatinine,  
SDMA, BUN

## Tubular function


Important in solute and water management

Urine concentration and protein; serum and urine electrolytes, glucose, acid-base

## Urine composition

Concentration/volume, pH, cellular and crystalline elements, infectious organisms, protein


## The functional nephron



# Renal tubules are where the action **really** is

## Tubular function

- + The actual work of the kidney primarily takes place here. Filtering, reabsorbing, and secreting solutes and water
- + Impact urine concentration and what is excreted
- + Dysfunction can impact electrolytes, protein levels, glucose, acid-base balance
- + Captured in chemistry panel and urinalysis





# Traditional renal **injury** markers are good, not great

- + Proteinuria
- + Hematuria, pyuria
- + Bacteriuria
- + Renal epithelial cells in the urine
- + Glucosuria (normoglycemia)
- + Cylindruria (casts)
- + Decreased USG

Granular casts



Cellular cast



Source: IDEXX SediVue Dx® images

Functional markers are  
in blood

Take-home message:

You can't assess kidney health  
without urine






Can we do any better?

YES!!! We can!

Cystatin B bridges the gap in  
our abilities to detect early  
and active renal injury





# What is cystatin B?




- + Member of cystatin family
  - + Protease inhibitors that help protect against leakage of proteolytic enzymes from lysosomes
  - + Trace amounts in the serum of healthy subjects
- + A small, intracellular protein
  - + 11 kDa (11,000 daltons)
  - + Ubiquitous in many cells, including proximal renal tubular cells
- + Freely filtered at the glomerulus
- + Increased urinary [cystatin B] indicates active, ongoing tubular injury
  - + Think of it as the ALT of the kidney


**Cystatin B** is a very small protein contained in epithelial cells of the renal tubules



During **active or acute kidney injury**, renal tubular epithelial cells can be damaged




# Kidney injury markers are additive to current indirect functional markers







# Kidney injury markers are additive to current indirect functional markers



# Traditional diagnostics only allow for case evaluation by functional markers



# Addition of an injury marker provides better case discrimination and management




Therefore, markers of tubular injury are earlier indicators of damage than functional markers

**By up to 2 days...**



# Cystatin B has value with evaluating patients with CKD as well!


- + CKD progressive and irreversible
- + Rate of progression unpredictable
- + Cystatin B identifies active, progressive injury in dogs with CKD
- + Increased urinary cystatin B (uCysB) in dogs with IRIS Stage 1 CKD predicts rapid progression
- + Identifies which dogs need more frequent monitoring



y-intercepts calculated from inverse urinary cystatin B (uCysB) vs. time



# Pathophysiology of AKI



## Common etiologies:

- + Infectious diseases
- + Nephrotoxins
- + Systemic diseases with secondary renal involvement (inflammation)
- + Alterations in hemodynamics
- + Obstructive disorders

# Causes of AKI include:

## Cat

- + Toxins (plants, chemotherapeutics)
- + Pyelonephritis
- + Acute pancreatitis
- + Marked dehydration
- + Obstructive disorders
- + Etiology unknown ~30%

## Dog

- + Toxins (plants, chemotherapeutics, foods)
- + Pyelonephritis
- + Acute pancreatitis
- + Marked dehydration
- + Obstructive disorders
- + Leptospirosis
- + Lyme nephritis
- + Congestive heart failure

# AKI can develop in hospitalized patients: Monitor and grade daily

- + Dehydration
- + Age > very young or old
- + Diuretic or nephrotoxic drug therapy
- + Hypokalemia or hypercalcemia
- + Sepsis
- + Congestive heart failure
- + Acute pancreatitis
- + Systemic hypertension
- + CKD

## **Avoid iatrogenic AKI!**

Nephrotoxic drugs  
Hemodynamic instability  
Fluid overload

## Initial treatment of AKI




# Initial therapy for AKI; back to the basics

- + Address hydration and volemic status
- + Institute disease-specific therapy whenever possible
- + Address complications of AKI
- + Avoid and monitor for complications of therapy
- + Keep the patient comfortable
- + Prevent further renal injury
- + Address nutrition

Use serial evaluations of IRIS AKI grading for objective measure of response to therapy; adapt/adjust accordingly

# Fluid therapy for kidney disease: Less may be more

- + Fluids are drugs—avoid overdose
- + Fluids do not improve kidney function
- + Hypervolemia causes AKI and kills patients that already have it
- + Not every patient with kidney disease (acute or chronic) needs fluids!!!





# Assessment of fluid therapy success is essential

- + Perfusion parameters: HR, CRT, mucous membranes, pulses, lactate, base excess
- + Body weight 2-4x/day: >5%–10% increase slow or stop fluids
- + Lung auscultation:  $\geq$  q12 hrs, more frequently if any changes in RR/RE



If azotemia worsens with IV fluid therapy, consider **decreasing** fluid rate.

Especially if total daily volume exceeds maintenance or if weight gain.





My dog ate some raisins.



Chemistry		3/14/24 3:32 AM		Collection		FREECATCH	
Glucose	105	63 - 114 mg/dL		Color	DARK YELLOW		
IDEXX SDMA	e 10	0 - 14 µg/dL		Clarity	TURBID		
Creatinine	1.0	0.5 - 1.5 mg/dL		Specific Gravity	1.049	>= 1.030	
IDEXX SDMA	10	0 - 14 µg/dL		pH	5.5	6.0 - 7.5	
Creatinine	1.0	0.5 - 1.5 mg/dL		Urine Protein	2+		
BUN	18	9 - 31 mg/dL		Glucose	NEGATIVE		
IDEXX Cystatin B (Urine)	>500	0 - 99 ng/mL		Ketones	NEGATIVE		
Potassium	5.1	4.0 - 5.4 mmol/L		Blood / Hemoglobin	3+		
Na: K Ratio	29	28 - 37		Bilirubin	1+		
Chloride	114	108 - 119 mmol/L		Urobilinogen	NORMAL		
TCO2 (Bicarbonate)	23	13 - 27 mmol/L		White Blood Cells	0-2		
Anion Gap	17	11 - 26 mmol/L		Red Blood Cells	10-15		
Total Protein	5.5	5.5 - 7.5 g/dL		Bacteria	RARE COCCI <9/HPF		
Albumin	3.0	2.7 - 3.9 g/dL		Additional Bacteria	RARE RODS <9/HPF		
				Casts	4+ (>10)/HPF		

# Rover

- + 3-year-old MC MixB
- + Confirmed raisin ingestion
- + Amount/time prior to presentation uncertain
- + Previously healthy

Elevated urinary cystatin B, proteinuria, cylindruria

# Three days later, after IV fluids for 48 hours

🧪

Chemistry

3/17/24

1:07 AM

📄

📖

IDEXX SDMA

8

0 - 14 µg/dL

📖

Creatinine

1.2

0.5 - 1.5 mg/dL

📖

BUN

26

9 - 31 mg/dL

📖

IDEXX Cystatin B (Urine)

<50

0 - 99 ng/mL

Ratio

📖

Phosphorus

5.1

2.5 - 6.1 mg/dL

📖

Calcium

9.5

8.4 - 11.8 mg/dL

📖

Sodium

148

142 - 152 mmol/L

📖

Potassium

5.1

4.0 - 5.4 mmol/L

📖

Na: K Ratio

29

28 - 37

📖

Chloride

114

108 - 119 mmol/L

📖

TCO2 (Bicarbonate)

25

13 - 27 mmol/L

📖

Anion Gap

14

11 - 26 mmol/L

📖

Total Protein

5.1

5.5 - 7.5 g/dL

📖

Albumin

2.7

2.7 - 3.9 g/dL

📖

Globulin

2.4

2.4 - 4.0 g/dL

📖

Glucose

NEGATIVE

📖

Ketones

a RACE

📖

Blood / Hemoglobin

3+

📖

Bilirubin

1+

📖

Urobilinogen

NORMAL

📖

White Blood Cells

0-2

📖

Red Blood Cells

30-50

📖

Bacteria

NONE SEEN

📖

Additional Bacteria

📖

Epithelial Cells

1+ (1-2)/HPF

📖

Mucus

NONE SEEN

📖

Casts

NONE SEEN

📖

Crystals

NONE SEEN

📖

Collection

FREECATCH

📖

J Vet Diagn Invest 17:223-231 (2005)

Canine renal pathology associated with grape or raisin ingestion: 10 cases

All dogs had degeneration or necrosis (or both) of proximal renal tubules with basement membranes remaining intact, and epithelial regeneration was observed in 5 out of 10 cases.

IDEXX

39

# Take-home

- + Acute kidney injury and chronic kidney disease are a continuum.
- + A COMPLETE urinalysis is of UTMOST importance when evaluating kidney (as well as systemic) disorders.
- + IDEALLY, patients at risk for renal injury (IRIS AKI grade I) are identified and managed BEFORE azotemia develops.
- + Fluid therapy paradigms have changed...dramatically.
- + Cystatin B, a *urine* biomarker, is a marker of ACTIVE renal tubular injury.







# **New Tests for Improved Early Diagnosis and Management of CKD**


Gregory F. Grauer, DVM, MS  
Diplomate, ACVIM (SAIM)  
Founding Member ACVNU  
Professor Emeritus  
Kansas State University

# CKD Diagnosis and Treatment:

## (Begin with the end in sight)


- ▶ Closer monitoring of SrCr and SDMA will facilitate early diagnosis of CKD
  - ▶ Longitudinal assessment of these parameters when available will provide better data than one-time evaluations
  - ▶ Monitoring serum SDMA concentrations will improve early diagnosis of CKD (2x more sensitive than SrCr)
- ▶ Use of FGF-23 will improve management of Phosphate imbalance/Phosphate overload
- ▶ Think “Inside the Reference Interval”

# **CKD is Irreversible and Often Progressive: No New Nephrons Can Be Produced**






# Potential effects of early detection and treatment of CKD



# Potential effects of early detection and treatment of CKD



# IRIS\* Classification of Canine and Feline CKD

<b>Creatinine SDMA**</b>	<b><u>Stage 1</u> Non-azotemic CKD</b>	<b><u>Stage 2</u> Normal to mild renal azotemia</b>	<b><u>Stage 3</u> Moderate renal azotemia</b>	<b><u>Stage 4</u> Severe renal azotemia</b>
<b>Creatinine Dogs/Cats (mg/dL)</b>	<b>&lt;1.4/&lt;1.6</b>	<b>1.4-2.8/1.6-2.8</b>	<b>2.9-5.0</b>	<b>&gt;5.0</b>
<b>SDMA Dogs/Cats (µg/dL)</b>	<b>&lt;18/&lt;18</b>	<b>18-35/18-25</b>	<b>36-54/26-38</b>	<b>&gt;54/&gt;38</b>
<b>Prevalence in Cats</b>	<b>33.3 %</b>	<b>37.2%</b>	<b>15.4 %</b>	<b>14.1 %</b>

Further sub-classify based on the presence or absence of proteinuria and systemic hypertension. \*[www.IRIS-kidney.com](http://www.IRIS-kidney.com)

*\*\*In the case of discrepancy between Creatinine and SDMA, consider body condition, and re-testing in 2-4 weeks. If values are persistent, consider assigning the patient to the higher stage.*



# IRIS\* Classification of Canine and Feline CKD

<b>Creatinine SDMA**</b>	<b><u>Stage 1</u> Non-azotemic CKD</b>	<b><u>Stage 2</u> Normal to mild renal azotemia</b>	<b><u>Stage 3</u> Moderate renal azotemia</b>	<b><u>Stage 4</u> Severe renal azotemia</b>
<b>Dogs/Cats (mg/dL)</b>	<1.4/<1.6	1.4-2.8/1.6-2.8	2.9-5.0	>5.0
<b>Dogs/Cats (µg/dL)</b>	<18/<18	18-35/18-25	36-54/26-38	>54/>38
<b>Prevalence in Cats</b>	33.3 %	37.2%	15.4 %	14.1 %

**Further sub-classify based on the presence or absence of proteinuria and systemic hypertension. \*[www.IRIS-kidney.com](http://www.IRIS-kidney.com)**

*\*\*In the case of discrepancy between Creatinine and SDMA, consider body condition, and re-testing in 2-4 weeks. If values are persistent, consider assigning the patient to the higher stage.*

# IRIS Guidelines: Diagnosis of Stage 1/Early Stage 2 CKD

- ▶ Kidney palpation or imaging abnormalities
- ▶ PU/PD due to loss of nephrons? (dogs > cats)
- ▶ Renal proteinuria (persistent w/ normal sediment)
- ▶ A persistent  $\uparrow$  in SDMA ( $\geq 14 \mu\text{g/dl}$ ) with sCr  $< 1.6 \text{ mg/dl}$  (cats) or  $< 1.4$  (dogs)
- ▶  $\uparrow$  in sCr within the RI without changes in muscle mass or hydration. Same for  $\uparrow$  in SDMA w/in RI
  - ▶ e.g., an increase in sCr from 0.7 to 1.4 mg/dl over several years could indicate  $\geq 50\%$  nephron loss.  $> 50\%$  because remaining nephrons have had time to undergo compensatory hypertrophy.



RIGHT K<sub>SU</sub>



# Serum Creatinine Concentration: Limitations


- ▶ Reflects muscle mass as well as GFR
- ▶ Influenced by methodology
  - ▶ Jaffe reaction vs. enzymatic; bench top vs. reference lab
- ▶ Variability in inter-laboratory reference ranges can lead to both false negative and false positive results for diagnosis of azotemia<sup>1,2</sup>
  - ▶ Longitudinal assessment (individualized baseline) with consistent methodology = best practice

1. Boozer L, et al. JVIM, 2002;16:354 (abstract)

2. Ulleberg T, et al. Acta Vet Scand 2011;53:25



# Relationship Between GFR and Serum Creatinine Concentration Is Not Linear






# Serum Symmetric Dimethylarginine (SDMA)

- ▶ Derived from intranuclear methylation of L-arginine by protein-arginine methyltransferases
- ▶ Released into circulation after proteolysis
- ▶ > 90% eliminated by glomerular filtration
  - ▶ Freely filtered with no tubular reabsorption
  - ▶ Non-renal influences appear to be minimal
- ▶ More sensitive than sCr: In two longitudinal studies, SDMA increased an average of 9 and 17 months prior to sCr, respectively in dogs and cats with CKD<sup>1,2</sup>


1. Hall JA, et al. JVIM 2014; 28:1676 (Feline)

2. Hall JA, et al. JVIM 2016; 30:794 (Canine)

# “Monkey” born 1996, FS, DSH




# “Monkey” born 1996, FS, DSH




CKD diagnosed 6/16/2010  
GFR 1.161 ml/min/kg  
40% loss from expected  
mean (1.94 ml/min/kg)

47% loss  
GFR 1.024 ml/min/kg

# “Monkey” born 1996, FS, DSH



# “Monkey” born 1996, FS, DSH



# Age-Specific Reference Intervals in Elderly Cats

Mortier F, et al JFMS 2023, Vol 25, Issue 11

- ▶ Using age-specific R is for sCr in mature adult (7-10 yrs) and senior cats ( > 10 yrs) improves health screening. (developed using ASVCP Guidelines)
- ▶ Standard RI for feline sCr was **0.9-2.3 mg/dl**
- ▶ Age-specific RI for mature adult cats would be **0.8-1.85 mg/dl**
- ▶ Age-specific RI for senior cats would be **0.7-1.86 mg/dl**



# Early Detection of CKD

- ▶ Serial determination of SrCr and SDMA concentrations
  - ▶ An increase  $\geq 0.3$  mg/dl is potentially real (vs. laboratory variation – when using the same lab/technique) for serum creatinine
  - ▶ An increase  $> 3.0$   $\mu$ g/dl is potentially real (vs. laboratory variation) for SDMA

# Summary of SCr and SDMA for Early Diagnosis of CKD


- ▶ SDMA is a more sensitive surrogate GFR marker than creatinine
- ▶ Longitudinal assessment of sCr and SDMA improves interpretation over single values
- ▶ Establishing baselines for subsequent longitudinal evaluation in individual patients is important
- ▶ Look for trends within the reference interval

# Renal Mineralization


Failure to control dietary phosphorus in dogs and cats with CKD can result in renal mineralization and fibrosis and progressive loss of nephrons

Ross, AJVR 43:1023, 1982

Finco, AJVR 53:2264, 1992








# IRIS Therapeutic Targets for Management of Hyperphosphatemia in CKD

IRIS Stage	Target Serum Phosphorus (mg/dl)	Management Options
1	2.5-4.5	Renal diet or Normal ration + binder
2	2.5-4.5	Renal diet $\pm$ binder
3	2.5-5.0	Renal diet + $\uparrow$ binder
4	2.5-6.0	Renal diet + $\uparrow\uparrow$ binder



# FGF-23 regulates Phos much like PTH regulates Ca



FGF-23, secreted from bone in response to  $\uparrow$  phosphorus. FGF-23  $\downarrow$  the density of renal Na-Phos cotransporters resulting in phosphaturia

# Utility of FGF-23 in Early-Stage CKD

- Plasma Phos  $>4.5$  mg/dl: Reduce Phos to  $<4.5$  for 30 days then measure FGF-23
- Plasma Phos  $<4.5$  mg/dl: Measure FGF-23
  - » FGF-23  $>400$  pg/ml: Restrict Phos intake and/or increase binder dose to  $\downarrow$  FGF-23
  - » FGF-23  $>300 <400$  pg/ml: Continue to monitor q 2-3 months
  - » FGF-23  $<300$  pg/ml: No need for further Phos restriction

# Utility of FGF-23 in Azotemic CKD

- Transition to a renal diet: Evaluate after 4-6 weeks
- Plasma Phos > Stage target range: Further reduce Phos intake and/or increase binder dose
  - » Phos within serum target range and FGF-23 >400 pg/ml: Further restrict Phos intake
  - » Phos within serum target range and FGF-23 <400 pg/ml: Continue current treatment

# Summary/Interpretation of Serum Phosphorus

- ▶ Phosphorus within the reference interval deserves more scrutiny – use IRIS Stage-specific Phos target guidelines
- ▶ Control with progressive use of renal therapeutic diets and enteric phosphate binders
- ▶ Once stage specific target intervals are achieved, further assess treatment efficacy with plasma FGF-23 concentrations





SS Phos  
Control



Don't just keep up.  
Stay ahead with  
*The Curated Mind*,  
IDEXX's monthly  
education newsletter.

# THE CURATED MIND

CONTINUING  
EDUCATION BROUGHT  
TO YOU BY IDEXX



LEARN EVOLVE REPEAT

**Scan the code** to stay in the know and receive updates about:

- + Upcoming live events
- + New on-demand education resources
- + The latest *Shake Up Your Workup* podcast
- + Articles from *The Vetiverse*
- + Upcoming in-person education opportunities





# THANK YOU

DON'T FORGET TO RATE YOUR SPEAKER  
AND SESSION IN THE APP!

Presented By

