Catalyst Dx* Chemistry Analyzer

QC/Maintenance Log

The Catalyst Dx*	analyzer performs a self-check every time it is powered on, wh	nen a
test is performed,	and when it is initialized.	

Year:			

Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Clean the internal components of the analyzer.												
Run VetTrol* quality control.												
Clean the centrifuge.												
Clean the waste drawer.												
Clean the fan filter.												
Clean the internal components of the analyzer.												
Clean the case.												
Clean the screen.												
Clean the waste drawer.												
	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.	Clean the internal components of the analyzer. Run VetTrol* quality control. Clean the centrifuge. Clean the waste drawer. Clean the fan filter. Clean the internal components of the analyzer. Clean the case. Clean the screen.

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





Catalyst One* Chemistry Analyzer

QC/Maintenance Log

The Cata	alyst One*	analyzer	performs	a self-checl	k every	time it is	powered	on,	when
a test is	performed	l, and whe	en it is init	ialized.					

Year:			

Frequency	Maintenance Required	January	February	March	April	Мау	June	July	August	September	October	November	December
Monthly	Clean the internal components of the analyzer.												
	Run VetTrol* quality control.												
	Clean the waste drawer.												
When Prompted	Clean the internal components of the analyzer.												
As Needed	Clean the case.												
	Clean the waste drawer.												

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





Coag Dx* Analyzer

QC/Maintenance Log

The Coag Dx* Analyzer performs a self-check every time it is activated, when a test is performed, and when an electronic quality control is run daily.

Month:							_																								
Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	:
As Needed	Clean the key																														

Required software upgrade (as provided)

pad, exterior case, and screen.

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





IDEXX Digital Cytology* Instrument

QC/Maintenance Log

Month:	

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Weekly	Shut down and restart the IDEXX Digital Cytology Instrument and its computer at least once per week.																															
Weekly	Change the stain and clean stain jars using a soft cloth and hot water (record the date the stain was changed on the coplin jar lids).																															

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





LaserCyte* Hematology Analyzer

QC/Maintenance Log

The LaserCyte analyzer is a liquid-based system that needs power
to maintain its cleansing cycle. If the LaserCyte analyzer is powered off for more
than 24 consecutive hours, please contact IDEXX Customer Support.

Quarter:	

Frequency	Maintenance Required	Week 1	Week 2	Week 3	Week 4	Week 5	Week 6	Week 7	Week 8	Week 9	Week 10	Week 11	Week 12
Weekly	Clean the air filter.												

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





ProCyte Dx* Hematology Analyzer

QC/Maintenance Log

Month:	

The analyzer will perform a background check when it restarts after exiting Standby mode. It takes approximately 8 minutes to exit Standby mode.

By default, the ProCyte Dx analyzer will set itself to Standby mode at 7:00 p.m. each day. To change the default setting, see the Operator's Guide or contact IDEXX Customer and Technical Support.

To power off the analyzer, tap **Turn Off ProCyte Dx** on the ProCyte Dx Instruments screen and follow the on-screen instructions. When powering on the IPU and ProCyte Dx analyzer, ensure that the IDEXX VetLab* Station is powered back on first.

Clean the fan filter while the analyzer is powered off for the weekly reboot.

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
As Needed	Change the reagent kit.																															
	Change the stain pack.																															
Weekly	Power off the analyzer and IPU.																															

Year:			

The Monthly Rinse procedure requires a 5% bleach solution of filtered, unscented bleach with distilled/deionized water. The ratio of bleach to water depends on the bleach concentration. (For example, Clorox* concentrated bleach has a 7.5% concentration, so the solution should be 2 parts bleach with 1 part deionized/distilled water.)

Use appropriate lab handling techniques when using e-CHECK* (XS) quality control.

Frequency	Maintenance Required	January	February	March	April	Мау	June	July	August	September	October	November	December
Monthly	Clean the fan filter.												
	Perform the Monthly Rinse procedure (when (prompted by the IDEXX VetLab Station).												
	Run e-CHECK (XS) L2 quality control.												

Required software upgrade (as provided)

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the ProCyte Dx analyzer, see the *ProCyte Dx Hematology Analyzer Operator's Guide* at idexx.com/product-info





ProCyte One* Hematology Analyzer

QC/Maintenance Log

Month:			

The ProCyte One* Hematology Analyzer performs most maintenance automatically outside of business hours. It runs a self-cleaning sequence after each sample run and cycles fluid daily outside of normal business hours. A more robust weekly maintenance sequence runs on Sundays at 3:00 a.m., including system checks, rinses, and an automated IDEXX SmartQC* control run. IDEXX will occasionally send software updates to the ProCyte One. No action is required on your part to install the updates.

IMPORTANT: The ProCyte One must remain powered on to receive software updates and to allow automated maintenance cycles to run. The analyzer will go into low-power mode after 15 minutes of idle time and will reboot periodically as needed during weekly maintenance.

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
As Needed	Change the reagent.																															
As Needed	Change the sheath.																															
As Needed	Change IDEXX SmartQC*.																															

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





IDEXX SediVue Dx* Urine Sediment Analyzer

QC/Maintenance Log

Month:	

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
As Needed	Empty the waste bin.																															
Weekly	Reboot the analyzer.								·				·																			

Year:			

Frequency	Maintenance Required	January	February	March	April	Мау	June	July	August	September	October	November	December
Monthly	Clean the analyzer.												
	Clean the fan filter.												
	Run L1 and L2 quality control.												

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





SNAP Pro* Analyzer

QC/Maintenance Log

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Weekly	Power the device off.																															
Weekly/As Needed	Clean the base, reflectors, and activator.																															

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





SNAPshot Dx* Analyzer

QC/Maintenance Log

Year:	

The SNAPshot Dx Analyzer calibrates automatically. It performs a self-check at the beginning of each test as well as an automatic systems check each time the analyzer is powered on.

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
As Needed	Clean the case.												
	Clean the SNAP* ports.												
	Clean the screen.												
	Clean the incubator wells.												

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





IDEXX VetAutoread* Hematology Analyzer

The IDEXX VetAutoread* analyzer performs a self-check every time it is powered on and whenever a Cal Rod is run.

QC/Maintenance Log

Month: _____

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily	Calibrate the analyzer with a Cal Rod.																															





IDEXX VetLab* Station

QC/Maintenance Log

Frequency	Maintenance Required	Week 1	Week 2	Week 3	Week 4	Week 5	Week 6	Week 7	Week 8	Week 9	Week 10	Week 11	Week 12
Weekly	Back up data.												
	Reboot the computer.												

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:





IDEXX VetLab* UA* Analyzer

QC/Maintenance Log

Month:									
--------	--	--	--	--	--	--	--	--	--

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily	At end of day, clean the test strip tray† with water and, if necessary, 70% isopropyl alcohol.																															
Weekly	Calibrate the analyzer every seven days.																										,					

†IMPORTANT: Be careful not to damage the gray reference pad during cleaning, and ensure that it is completely clean and dry before reading a test strip.





VetStat* Electrolyte and Blood Gas Analyzer

The VetStat analyzer performs a self-check every time it is activated, when a test is performed, and when an electronic quality control is run daily.

QC/Maintenance Log

Month:	

Frequency	Maintenance Required	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily	Run a standard reference cassette (SRC), alternate levels.																															
Weekly	Clean the sample measurement chamber.																															

Year:		

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Run OPTI* CHECK quality control (alternate levels).												
Quarterly	Run a hemoglobin calibration cassette (HbCC) and reinstall the standard reference cassettes (i.e., rescan SRC setup barcodes).												
Annually	Replace peristaltic pump cartridge.												
	Change gas I/O port.												





VetTest* Chemistry Analyzer

QC/Maintenance Log

Year:						

Frequency	Maintenance Required	January	February	March	April	May	June	July	August	September	October	November	December
Monthly	Run VetTrol* quality control with the quality control panel.												
	Clean pipettor actuator collar with												

Required software upgrade (as provided)

an alcohol swab.

New version number:	New version number:	New version number:	New version number:
Date uploaded:	Date uploaded:	Date uploaded:	Date uploaded:

For more information about maintaining the VetTest Analyzer, see the VetTest Chemistry Analyzer Operator's Guide at idexx.com/product-info.





The VetTest analyzer performs a self-check every time it is activated, when a test is

performed, and when an electronic quality control is run daily.