


Catalyst SDMA Test

Quick Reference Guide

Storage and handling requirements

- + Store in the freezer (do not store in the freezer door).
- + No warming required—run directly from the freezer.
- + Catalyst™ SDMA tests can be stored at room temperature for up to 8 hours. After 8 hours, store SDMA tests in the freezer.
- + Sample volume required:
 - Whole blood: 0.7 cc (700 µL)
 - Serum/plasma: 100 µL (65 µL minimum; 300 µL if running with other slides)

Just load and go!



Frequently asked questions

Question	Answer						
Do I need to run anything with the SDMA slide?	SDMA can be run alone as a single slide.						
What sample types can be run on the SDMA test?	Compatible sample types include serum, lithium heparin plasma, and whole blood using the Catalyst™ Lithium Heparin Whole Blood Separator.						
Can I run previously frozen samples on the SDMA test?	Yes, previously frozen samples can be used when running the Catalyst SDMA Test.						
Which species have been validated for the SDMA test and what are the numerical results?	<table border="1"> <thead> <tr> <th>Species</th> <th>U.S. Units/ S.I./French S.I. Units</th> </tr> </thead> <tbody> <tr> <td>Canine</td> <td>0–100 µg/dL</td> </tr> <tr> <td>Feline</td> <td>0–100 µg/dL</td> </tr> </tbody> </table>	Species	U.S. Units/ S.I./French S.I. Units	Canine	0–100 µg/dL	Feline	0–100 µg/dL
Species	U.S. Units/ S.I./French S.I. Units						
Canine	0–100 µg/dL						
Feline	0–100 µg/dL						
Can SDMA be run with other slides?	<p>Yes! The test can be run alone or with other slides as part of a comprehensive patient profile. For example, you could run a Catalyst™ Lyte 4 CLIP, a Catalyst™ Chem 17 CLIP, a Catalyst™ Total T₄ Test, and a Catalyst™ SDMA Test with one patient sample.</p> <p>IMPORTANT: The SDMA slide should not be run with a Catalyst™ Chem 18, Chem 16, Chem 11, or NSAID7 CLIP (an SDMA slide is already included in these CLIPs).</p>						
When running with other slides, should the SDMA test be loaded in a particular order?	<ul style="list-style-type: none"> + The SDMA slide can be run in any order after electrolytes. + When running SDMA with electrolyte slides, always load the electrolyte slides first. + For the quickest time to results for SDMA, the recommended load order is Lyte 4 CLIP on the bottom, followed by a chemistry CLIP (e.g., Chem 17, Chem 10, etc.), SDMA, any additional slides, and TT₄ on top. 						
What is the run time for the SDMA test?	SDMA results are available ~7 minutes after the start of the run.						
Is it okay to dilute samples that are undergoing SDMA testing?	IDEXX does not support the use of diluted samples when running the Catalyst SDMA Test.						
How often can SDMA tests be left at room temperature and then returned to the freezer?	Once at room temperature, SDMA tests can be returned to the freezer up to 5 times as long as they are unopened.						
What if an SDMA test is accidentally refrigerated?	Any test suspected of being refrigerated should be discarded.						