



JANUARY 17-21 • NAVC.COM • ORLANDO, FL

WORLD CLASSIC

CELEBRATING THE CHAMPIONS OF CARE

VMX
VETERINARY MEETING & EXPO



CBC edu'CAT'ion: Feline hematology Secrets

Fred L. Metzger DVM, MRCVS, DABVP

Dennis B. DeNicola DVM, PhD, DACVP

Alan H. Rebar DVM, PhD, DACVP




Old dogs, old tricks

- Feline Reticulocytes
- Feline Heinz Bodies
- Identification of band/immature neutrophils-blood film
- Cytograms/dot plots


Feline reticulocytes

- 0.4% present in normal health (aggregate)
- Two types of reticulocytes:
 - Aggregate:
 - ↓ ~12–24 hours
 - Punctate:
 - ↓ ~ 10–14 days
- Mature RBC



Definition of regeneration in cats


- How many reticulocytes?
- Only aggregate reticulocytes?
- Aggregate and punctate reticulocytes?



Sources:


1. Cramer DV, Lewis RM. Reticulocyte response in the cat. *JAVMA*, 1972;160(1):61-67.
2. Alsaker RD, Laber J, Stevens JB, Perman V. A comparison of polychromasia and reticulocyte counts in assessing erythrocytic regenerative response in the cat. *JAVMA*. 1977;170(1):39-41.
3. Fan LC, Dorner JL, Hoffman WE. Reticulocyte response and maturation in experimental acute blood loss anemia in the cat. *J Am Anim Hosp Assoc*. 1978;14:219-224.

Feline *Mycoplasma haemofelis* anemia



Hospital day 2 PCV = 8%

Feline *Mycoplasma haemofelis* anemia



Hospital day 10 PCV = 22%

Heinz Bodies

- **In cats, Heinz bodies are particularly significant for a few reasons:**
- **Cats are more prone to Heinz body formation** than many other species. This is because feline hemoglobin has more sulfhydryl groups, making it more susceptible to oxidative damage.
- **Small numbers of Heinz bodies (up to 5%) can be normal** in healthy cats and may not indicate disease.
- **Increased numbers of Heinz bodies** can be seen with exposure to oxidative agents, such as certain drugs (e.g., acetaminophen), foods (e.g., onions, garlic), or diseases (e.g., diabetes mellitus, lymphoma, hyperthyroidism).
- **Heinz bodies can cause hemolytic anemia** if present in large numbers, as affected red blood cells are more likely to be removed by the spleen.

Poikilocytes in cats


- Spherocytes
- Acanthocytes

Inflammation


- Dohle bodies
 - Döhle bodies in cats as in other species are cytoplasmic inclusions found in neutrophils, appearing as pale blue, angular structures
 - Often associated with toxic change during inflammation in most species
 - In cats – often identified in neutrophils from healthy cats
 - May be seen during inflammation but difficult to differentiate
 - Need to look for other indicators of toxicity (cytoplasmic basophilia, foamy to vacuolated cytoplasm, etc.

Normal Canine and Feline WBC Dot Plots

Normal WBC Dot Plot (Canine)



Normal WBC Dot Plot (Feline)



New Tricks from Old Dogs

- Reticulocytosis without anemia
- Reticulocyte Hemoglobin
- Feline SAA- Serum Amyloid A
- Red Cell and Platelet Indices
- AI Morphology
- Feline Pancreatic Lipase

Feline acute phase proteins

- Major APP in the cat
 - Serum amyloid A (SAA)
- SAA sample
 - Serum
 - Heparinized and EDTA plasma
 - Body cavity fluids (synovial, peritoneal)
- SAA stability
 - Room temp preferred over refrigerated
 - >24-hour storage freeze sample
 - Frozen 3 months

| Species | Major APP | Moderate APP | Minor APP |
|---------|-----------|--------------|-----------|
| Feline | SAA | AGP | Hp |
| | | CRP | Cp |


SAA= Serum Amyloid A; AGP= alpha-1 acid glycoprotein; Hp= Haptoglobin;
CRP= C Reactive Protein; Cp= ceruloplasmin

Platelets by the numbers

- Thrombogram
 - PLT / PCT – count / mass

| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------|----------------|--------------------|-----|--------|------|
| PLT | 225 K/ μ L | 148 - 484 | | | |
| MPV | 12.6 fL | 8.7 - 13.2 | | | |
| PDW | 12.3 fL | 9.1 - 19.4 | | | |
| PCT | 0.28 % | 0.14 - 0.46 | | | |

Hematocrit tube




Analogous to HCT

$$\text{HCT (\%)} = (\text{RBC} \times \text{MCV}) / 10$$

Platelets by the numbers


- Thrombogram
 - PLT / PCT – count / mass

| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------|----------------|--------------------|-----|--------|------|
| PLT | 225 K/ μ L | 148 - 484 | | | |
| MPV | 12.6 fL | 8.7 - 13.2 | | | |
| PDW | 12.3 fL | 9.1 - 19.4 | | | |
| PCT | 0.28 % | 0.14 - 0.46 | | | |





Platelets by the numbers

- Thrombogram
 - PLT / PCT – count / mass
 - MPV – platelet volume




| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------|----------------|--------------------|-----|--------|------|
| PLT | 225 K/ μ L | 148 - 484 | | | |
| MPV | 12.6 fL | 8.7 - 13.2 | | | |
| PDW | 12.3 fL | 9.1 - 19.4 | | | |
| PCT | 0.28 % | 0.14 - 0.46 | | | |

Species variation in platelet size



Platelets by the numbers

- Thrombogram
 - PLT / PCT – count / mass
 - MPV – platelet volume
 - PDW – variation in platelet size



| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------|----------------|--------------------|-----|--------|------|
| PLT | 225 K/ μ L | 148 - 484 | | | |
| MPV | 12.6 fL | 8.7 - 13.2 | | | |
| PDW | 12.3 fL | 9.1 - 19.4 | | | |
| PCT | 0.28 % | 0.14 - 0.46 | | | |

Clinical Case Study

- **History/Presentation**

- Adopted 5 months earlier, soon to receive final vaccinations
- In-door cat since adoption
- Presented with sudden onset lethargy, anorexia and ataxia

- **Physical examination**

- Temp: 102.1° F, Pulse: 220 bpm, Resp: 80
- Pale pink mucous membranes, CRT < 2 sec
- Dull mentation / Depressed
- Walking slowly with ataxia
- Eventual blindness

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

- Initial laboratory evaluation three weeks before present CBC
 - Mild decreased MCV – 33.3 fL (RI; 35.9 – 53.1)
 - Marked decreased PLT – 13 K/ μ L (RI; 151 – 600)
 - Mild ALP increase – 154 U/L (RI; 14 – 111)
 - Mild PHOS increase – 7.8 mg/dL (RI; 3.1 – 7.5)
 - Marked Ammonia increase – 217 μ mol/L (RI; 0 – 95)
- Additional diagnostics
 - Ultrasound / CT scan of abdomen – possible extrahepatic shunt
- Patient management
 - Surgical correction of shunt

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|-------------------------------------|------------------|--------------------|------|--------|------|
| ProCyte Dx (March 12, 2024 8:19 AM) | | | | | |
| RBC | 9.91 M/ μ L | 6.54 - 12.20 | | | |
| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
| ProCyte Dx (March 12, 2024 8:19 AM) | | | | | |
| PLT | 29 K/ μ L | 151 - 600 | LOW | | |
| MPV | 14.7 fL | 11.4 - 21.6 | | | |
| PCT | 0.04 % | 0.17 - 0.86 | LOW | | |
| WBC | 14.45 K/ μ L | 2.87 - 17.02 | | | |
| %NEU | 48.5 % | | | | |
| %LYM | 38.6 % | | | | |
| %MONO | 2.5 % | | | | |
| %EOS | 0.6 % | | | | |
| %BASO | 9.8 % | | | | |
| NEU | 7.02 K/ μ L | 2.30 - 10.29 | | | |
| LYM | 5.58 K/ μ L | 0.92 - 6.88 | | | |
| MONO | 0.36 K/ μ L | 0.05 - 0.67 | | | |
| EOS | 0.08 K/ μ L | 0.17 - 1.57 | LOW | | |
| BASO | 1.41 K/ μ L | 0.01 - 0.26 | HIGH | | |
| PLT | 29 K/ μ L | 151 - 600 | LOW | | |
| MPV | 14.7 fL | 11.4 - 21.6 | | | |
| PCT | 0.04 % | 0.17 - 0.86 | LOW | | |

ProCyte Dx is a trademark or registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH




Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|-------------------------------------|---------------|--------------------|-----|--------|------|
| ProCyte Dx (March 12, 2024 8:19 AM) | | | | | |
| PLT | 29 K/ μ L | 151 - 600 | LOW | | |
| MPV | 14.7 fL | 11.4 - 21.6 | | | |
| PCT | 0.04 % | 0.17 - 0.86 | LOW | | |



Decreased PLT
Decreased PCT



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat


Clinical Case Study

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------------------------------------|---------------|--------------------|-----|--------|------|
| ProCyt Dx (March 12, 2024 8:19 AM) | | | | | |
| PLT | 29 K/ μ L | 151 - 600 | LOW | | |
| MPV | 14.7 fL | 11.4 - 21.6 | | | |
| PCT | 0.04 % | 0.17 - 0.86 | LOW | | |

How can we confirm the thrombocytopenia?




Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study


| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------|---------|--------------------|-----|--------|------|
|------|---------|--------------------|-----|--------|------|

ProCyt Dx (March 12, 2024 8:19 AM)

| | | | | | |
|-----|---------------|-------------|-----|--|--|
| PLT | 29 K/ μ L | 151 - 600 | LOW | | |
| MPV | 14.7 fL | 11.4 - 21.6 | | | |
| PCT | 0.04 % | 0.17 - 0.86 | LOW | | |



■ RBC ■ RETICS ■ PLT ■ RBC Frags ■ WBC



■ RBC ■ RETICS ■ PLT ■ RBC Frags ■ WBC

- Overall identification of different cell populations appears good – you can trust the data
- Platelet digitized events are low compared to normal – supports thrombocytopenia
Continuum of digitized events suggests very large or clumps of platelets

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

Boogeyman


- 7 months old
- Neutered male
- Black and White DSH

| Test | Results | Reference Interval | LOW | NORMAL | HIGH |
|------|---------|--------------------|-----|--------|------|
|------|---------|--------------------|-----|--------|------|


ProCyte Dx (March 12, 2024 8:19 AM)

| | | | | | |
|-----|---------------|-------------|-----|--|--|
| PLT | 29 K/ μ L | 151 - 600 | LOW | | |
| MPV | 14.7 fL | 11.4 - 21.6 | | | |
| PCT | 0.04 % | 0.17 - 0.86 | LOW | | |

Strongly supports platelet clumping



■ NEU ■ LYM ■ MONO ■ EOS ■ BASO ■ URBC



■ NEU ■ LYM ■ MONO ■ EOS ■ BASO ■ URBC




- Some significant variation from normal pattern – needs blood film review
 - Extremely low numbers of eosinophils
 - Unlysed RBCs and NEUs melding together
- Platelet clumping oftentimes presents as a curvilinear array of digitized events extending from the lower left corner upward to the right



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

Examples of PLT clumping with the ProCyt Dx



Boogeyman





- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

Example of PLT clumping with the ProCyt One



Boogeyman


- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

Feathered edge of blood film



Boogeyman


- 7 months old
- Neutered male
- Black and White DSH



Clinical Case Study

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH




Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

Clinical Case Study

Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

Monolayer of blood film



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

IDEXX inVue Dx™ Cellular Analyzer

Revolutionary workflow

Slide-free, load-and-go, results in 10 minutes

Plug-and-play integration with VetConnect® PLUS and IDEXX VetLab® Station

Give valuable time back for patient care

Deeper insights

Multidimensional interrogation of cells in natural state


Multiple wavelengths of light and fluorescence

Board-certified pathologist trained





IDEXX inVue Dx, VetConnect PLUS, and IDEXX VetLab Station are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

IDEXX inVue Dx™ Cellular Analyzer




Platelet Clumping



$$\text{PCT (\%)} = (\text{PLT} \times \text{MPV}) / 10,000$$

Approach to the evaluation of platelets



The Laboratory Retriever Team



Alan H. Rebar, DVM, PhD, DACVP



Dennis B. DeNicola, DVM, PhD, DACVP



Dr. Fred Metzger, DVM, MRCVS, DABVP

Catalyst New Menu Consistently Drives Adoption at Scale

% of **U.S. Catalyst** customers **adopted in the first year of launch***

>65%



Catalyst® Lyte 4 CLIP

Launched December 2008

>55%



Catalyst® Total T4 Test

Launched February 2015

>50%



Catalyst® SDMA

Launched January 2018

>50%



Catalyst®
Pancreatic Lipase

Launched September 2024

*Customers who adopted in first year of launch had run at least one consumable slide on Catalyst Dx® Chemistry Analyzer or Catalyst One® Chemistry Analyzer in the last 30 days of 12 months post-launch for Lyte 4 CLIP, Total T4 and SDMA, and at least one

consumable slide on Catalyst Dx® Chemistry Analyzer or Catalyst One® Chemistry Analyzer in the last 30 days as of 7/30/25 for Pancreatic Lipase.

© 2025 IDEXX Laboratories, Inc. All rights reserved.

Don't just keep up.
Stay ahead with
The Curated Mind,
IDEXX's monthly
education newsletter.

THE CURATED MIND

CONTINUING
EDUCATION BROUGHT
TO YOU BY IDEXX



LEARN EVOLVE REPEAT

Scan the code to stay in the know and receive updates about:

- + Upcoming live events
- + New on-demand education resources
- + The latest *Shake Up Your Workup* podcast
- + Articles from *The Vetiverse*
- + Upcoming in-person education opportunities





THANK YOU

DON'T FORGET TO RATE YOUR SPEAKER
AND SESSION IN THE APP!

Presented By

