

# CBC edu'CAT'ion: Feline hematology Secrets

Fred L. Metzger DVM, MRCVS, DABVP

Dennis B. DeNicola DVM, PhD, DACVP

Alan H. Rebar DVM, PhD, DACVP




# Old dogs, old tricks

- Feline Reticulocytes
- Feline Heinz Bodies
- Identification of band/immature neutrophils-blood film
- Cytograms/dot plots


# Feline reticulocytes

- 0.4% present in normal health (aggregate)
- Two types of reticulocytes:
  - Aggregate:
    - ↓ ~12–24 hours
  - Punctate:
    - ↓ ~ 10–14 days
  - Mature RBC



# Definition of regeneration in cats


- How many reticulocytes?
- Only aggregate reticulocytes?
- Aggregate and punctate reticulocytes?



#### Sources:


1. Cramer DV, Lewis RM. Reticulocyte response in the cat. *JAVMA*, 1972;160(1):61-67.
2. Alsaker RD, Laber J, Stevens JB, Perman V. A comparison of polychromasia and reticulocyte counts in assessing erythrocytic regenerative response in the cat. *JAVMA*. 1977;170(1):39-41.
3. Fan LC, Dorner JL, Hoffman WE. Reticulocyte response and maturation in experimental acute blood loss anemia in the cat. *J Am Anim Hosp Assoc*. 1978;14:219-224.

# Feline *Mycoplasma haemofelis* anemia



Hospital day 2    PCV = 8%

# Feline *Mycoplasma haemofelis* anemia



Hospital day 10 PCV = 22%

# Heinz Bodies

- **In cats, Heinz bodies are particularly significant for a few reasons:**
- **Cats are more prone to Heinz body formation** than many other species. This is because feline hemoglobin has more sulfhydryl groups, making it more susceptible to oxidative damage.
- **Small numbers of Heinz bodies (up to 5%) can be normal** in healthy cats and may not indicate disease.
- **Increased numbers of Heinz bodies** can be seen with exposure to oxidative agents, such as certain drugs (e.g., acetaminophen), foods (e.g., onions, garlic), or diseases (e.g., diabetes mellitus, lymphoma, hyperthyroidism).
- **Heinz bodies can cause hemolytic anemia** if present in large numbers, as affected red blood cells are more likely to be removed by the spleen.

# Poikilocytes in cats


- Spherocytes
- Acanthocytes

# Inflammation


- Dohle bodies
  - Döhle bodies in cats as in other species are cytoplasmic inclusions found in neutrophils, appearing as pale blue, angular structures
  - Often associated with toxic change during inflammation in most species
  - In cats – often identified in neutrophils from healthy cats
    - May be seen during inflammation but difficult to differentiate
    - Need to look for other indicators of toxicity (cytoplasmic basophilia, foamy to vacuolated cytoplasm, etc.)

# Normal Canine and Feline WBC Dot Plots

Normal WBC Dot Plot (Canine)



Normal WBC Dot Plot (Feline)



# New Tricks from Old Dogs

- Reticulocytosis without anemia
- Reticulocyte Hemoglobin
- Feline SAA- Serum Amyloid A
- Red Cell and Platelet Indices
- AI Morphology
- Feline Pancreatic Lipase

# Feline acute phase proteins

- Major APP in the cat
  - Serum amyloid A (SAA)
- SAA sample
  - Serum
  - Heparinized and EDTA plasma
  - Body cavity fluids (synovial, peritoneal)
- SAA stability
  - Room temp preferred over refrigerated
    - >24-hour storage freeze sample
  - Frozen 3 months

Species	Major APP	Moderate APP	Minor APP
Feline	SAA	AGP	Hp
		CRP	Cp


SAA= Serum Amyloid A; AGP= alpha-1 acid glycoprotein; Hp= Haptoglobin;  
CRP= C Reactive Protein; Cp= ceruloplasmin

# Platelets by the numbers

- Thrombogram
  - PLT / PCT – count / mass

Test	Results	Reference Interval	LOW	NORMAL	HIGH
PLT	225 K/ $\mu$ L	148 - 484		█	
MPV	12.6 fL	8.7 - 13.2		█	
PDW	12.3 fL	9.1 - 19.4		█	
PCT	0.28 %	0.14 - 0.46		█	

Hematocrit tube




Analogous to HCT

$$\text{HCT (\%)} = (\text{RBC} \times \text{MCV}) / 10$$

# Platelets by the numbers


- Thrombogram
  - PLT / PCT – count / mass

Test	Results	Reference Interval	LOW	NORMAL	HIGH
PLT	225 K/ $\mu$ L	148 - 484		█	
MPV	12.6 fL	8.7 - 13.2			█
PDW	12.3 fL	9.1 - 19.4		█	
PCT	0.28 %	0.14 - 0.46		█	







# Platelets by the numbers

- Thrombogram
  - PLT / PCT – count / mass
  - MPV – platelet volume




Test	Results	Reference Interval	LOW	NORMAL	HIGH
PLT	225 K/ $\mu$ L	148 - 484			
MPV	12.6 fL	8.7 - 13.2			
PDW	12.3 fL	9.1 - 19.4			
PCT	0.28 %	0.14 - 0.46			

# Species variation in platelet size



# Platelets by the numbers

- Thrombogram
  - PLT / PCT – count / mass
  - MPV – platelet volume
  - PDW – variation in platelet size



Test	Results	Reference Interval	LOW	NORMAL	HIGH
PLT	225 K/ $\mu$ L	148 - 484			
MPV	12.6 fL	8.7 - 13.2			
PDW	12.3 fL	9.1 - 19.4			
PCT	0.28 %	0.14 - 0.46			

# Clinical Case Study

- **History/Presentation**

- Adopted 5 months earlier, soon to receive final vaccinations
- In-door cat since adoption
- Presented with sudden onset lethargy, anorexia and ataxia

- **Physical examination**

- Temp: 102.1° F, Pulse: 220 bpm, Resp: 80
- Pale pink mucous membranes, CRT < 2 sec
- Dull mentation / Depressed
- Walking slowly with ataxia
- Eventual blindness

## **Boogeyman**

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study

- Initial laboratory evaluation three weeks before present CBC
  - Mild decreased MCV – 33.3 fL (RI; 35.9 – 53.1)
  - Marked decreased PLT – 13 K/ $\mu$ L (RI; 151 – 600)
  - Mild ALP increase – 154 U/L (RI; 14 – 111)
  - Mild PHOS increase – 7.8 mg/dL (RI; 3.1 – 7.5)
  - Marked Ammonia increase – 217  $\mu$ mol/L (RI; 0 – 95)
- Additional diagnostics
  - Ultrasound / CT scan of abdomen – possible extrahepatic shunt
- Patient management
  - Surgical correction of shunt

## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study

Test	Results	Reference Interval	LOW	NORMAL	HIGH
ProCyte Dx (March 12, 2024 8:19 AM)					

RBC	9.91 M/ $\mu$ L	6.54 - 12.20			
-----	-----------------	--------------	--	--	--

Test	Results	Reference Interval	LOW	NORMAL	HIGH
ProCyte Dx (March 12, 2024 8:19 AM)					

PLT	29 K/ $\mu$ L	151 - 600	LOW		
MPV	14.7 fL	11.4 - 21.6			
PCT	0.04 %	0.17 - 0.86	LOW		

WBC	14.45 K/ $\mu$ L	2.87 - 17.02			
-----	------------------	--------------	--	--	--

%NEU	48.5 %				
%LYM	38.6 %				
%MONO	2.5 %				
%EOS	0.6 %				
%BASO	9.8 %				

NEU	7.02 K/ $\mu$ L	2.30 - 10.29			
LYM	5.58 K/ $\mu$ L	0.92 - 6.88			
MONO	0.36 K/ $\mu$ L	0.05 - 0.67			
EOS	0.08 K/ $\mu$ L	0.17 - 1.57	LOW		
BASO	1.41 K/ $\mu$ L	0.01 - 0.26	HIGH		
PLT	29 K/ $\mu$ L	151 - 600	LOW		
MPV	14.7 fL	11.4 - 21.6			
PCT	0.04 %	0.17 - 0.86	LOW		

ProCyte Dx is a trademark or registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study


## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

Test	Results	Reference Interval	LOW	NORMAL	HIGH
------	---------	--------------------	-----	--------	------

ProCyt Dx (March 12, 2024 8:19 AM)

PLT	29 K/ $\mu$ L	151 - 600	LOW		
MPV	14.7 fL	11.4 - 21.6			
PCT	0.04 %	0.17 - 0.86	LOW		



Decreased PLT  
Decreased PCT



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat


# Clinical Case Study

## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

Test	Results	Reference Interval	LOW	NORMAL	HIGH
ProCyte Dx (March 12, 2024 8:19 AM)					
PLT	29 K/ $\mu$ L	151 - 600	LOW		
MPV	14.7 fL	11.4 - 21.6			
PCT	0.04 %	0.17 - 0.86	LOW		

How can we confirm the thrombocytopenia?



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study



## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

Test	Results	Reference Interval	LOW	NORMAL	HIGH
------	---------	--------------------	-----	--------	------

ProCyte Dx (March 12, 2024 8:19 AM)

PLT	29 K/ $\mu$ L	151 - 600	LOW		
MPV	14.7 fL	11.4 - 21.6			
PCT	0.04 %	0.17 - 0.86	LOW		



- Overall identification of different cell populations appears good – you can trust the data
- Platelet digitized events are low compared to normal – supports thrombocytopenia
- Continuum of digitized events suggests very large or clumps of platelets



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study

## Boogeyman



- 7 months old
- Neutered male
- Black and White DSH

Test	Results	Reference Interval	LOW	NORMAL	HIGH
------	---------	--------------------	-----	--------	------

ProCyte Dx (March 12, 2024 8:19 AM)

PLT	29 K/ $\mu$ L	151 - 600	LOW		
MPV	14.7 fL	11.4 - 21.6			
PCT	0.04 %	0.17 - 0.86	LOW		

Strongly supports platelet clumping






- Some significant variation from normal pattern – needs blood film review
  - Extremely low numbers of eosinophils
  - Unlysed RBCs and NEUs melding together
- Platelet clumping oftentimes presents as a curvilinear array of digitized events extending from the lower left corner upward to the right



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study

Examples of PLT clumping with the ProCyte Dx



## Boogeyman





- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study

Example of PLT clumping with the ProCyt One




## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study



## Boogeyman


- 7 months old
- Neutered male
- Black and White DSH



# Clinical Case Study

## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH




Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# Clinical Case Study

## Boogeyman

- 7 months old
- Neutered male
- Black and White DSH

Monolayer of blood film



Not Boogeyman: only a representation of a black and white Domestic Shorthair cat

# IDEXX inVue Dx™ Cellular Analyzer

## Revolutionary workflow

---

Slide-free, load-and-go, results in 10 minutes

Plug-and-play integration with VetConnect® PLUS and IDEXX VetLab® Station

Give valuable time back for patient care

## Deeper insights

---

Multidimensional interrogation of cells in natural state


Multiple wavelengths of light and fluorescence

Board-certified pathologist trained





IDEXX inVue Dx, VetConnect PLUS, and IDEXX VetLab Station are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

# IDEXX inVue Dx™ Cellular Analyzer




Platelet Clumping



$$\text{PCT (\%)} = (\text{PLT} \times \text{MPV}) / 10,000$$



# Approach to the evaluation of platelets



# The Laboratory Retriever Team



Alan H. Rebar, DVM, PhD, DACVP



Dennis B. DeNicola, DVM, PhD, DACVP



Dr. Fred Metzger, DVM, MRCVS, DABVP

# Catalyst New Menu Consistently Drives Adoption at Scale

% of **U.S. Catalyst** customers **adopted in the first year of launch\***

>65%



Catalyst® Lyte 4 CLIP

*Launched December 2008*

>55%



Catalyst® Total T4 Test

*Launched February 2015*

>50%



Catalyst® SDMA

*Launched January 2018*

>50%



Catalyst®  
Pancreatic Lipase

*Launched September 2024*

\*Customers who adopted in first year of launch had run at least one consumable slide on Catalyst Dx® Chemistry Analyzer or Catalyst One® Chemistry Analyzer in the last 30 days of 12 months post-launch for Lyte 4 CLIP, Total T4 and SDMA, and at least one

consumable slide on Catalyst Dx® Chemistry Analyzer or Catalyst One® Chemistry Analyzer in the last 30 days as of 7/30/25 for Pancreatic Lipase.

© 2025 IDEXX Laboratories, Inc. All rights reserved.

Don't just keep up.  
Stay ahead with  
*The Curated Mind*,  
IDEXX's monthly  
education newsletter.

# THE CURATED MIND

CONTINUING  
EDUCATION BROUGHT  
TO YOU BY IDEXX



LEARN EVOLVE REPEAT

Scan the code to stay in the know and receive updates about:

- + Upcoming live events
- + New on-demand education resources
- + The latest *Shake Up Your Workup* podcast
- + Articles from *The Vetiverse*
- + Upcoming in-person education opportunities

