

Digital Imaging, No Delays

All it takes is a little preparation to ensure a smooth installation of your new IDEXX I-Vision Digital Imaging system



Congratulations—upgrading your practice with the IDEXX I-Vision Digital Imaging system will help you deliver the highest standard of patient care. This quick guide will help you avoid common pitfalls and get your system up and running quickly.

#1: Be sure there are no technical issues with your practice generator or electrical system.

Prior to installation, it's important that your x-ray generator is ready to use and has been calibrated by your local service technician within the last 24 months.

Your new IDEXX I-Vision Digital Imaging System requires a standard, properly grounded 3-pronged power outlet.

#2: Clear adequate space for your new digital system.

In addition to removing your old processor and equipment, you'll need enough space and furniture to support the IDEXX I-Vision CR system scanner (if applicable) and the IDEXX-PACS capture station computer, monitor, keyboard and mouse.

If you don't have enough space, a cart will be needed. Please make sure the cart is assembled prior to your digital installation.

#3: Prepare your practice network.

An active network connection must be installed, tested (CAT5 or CAT6 cable) and connected to the server within 5 feet of the capture station.

An Internet connection will ensure you can fully utilize all system features and support.

#4: Upgrade your practice information management software (PIMS).

Cornerstone® software users should upgrade to the latest version prior to digital installation. You will need to load a new activation key, which will be provided after the preinstallation kick-off with IDEXX.

Other practice management systems may also require an upgrade to maximize the IDEXX I-Vision Digital Imaging system's integration capabilities. Contact your PIMS provider for information on how to upgrade prior to digital installation.

#5: Make sure your staff is available for training.

Learning the ins and outs of your new IDEXX I-Vision system is easy with the right training. To make the most of your system, be sure that at least two of your staff members are available to be trained by the on-site installer on the day(s) of your digital installation.

Questions?

Contact your IDEXX representative or call 1-877-433-9948.