



Establish a protocol to improve reproductive efficiency

Collaborate to identify open animals early.

Producers, veterinarians, and laboratories that work together can act early to reduce open days and address health issues that impact pregnancy.

Trust the science of specific detection.

Pregnancy-associated glycoproteins (PAGs) are only produced in the presence of an embryo or fetus.

Depend on accurate results.

IDEXX pregnancy tests are highly specific for PAGs, providing definitive results to supplement current methods and guide veterinary recommendations for animal health and productivity.

.....

Contact your IDEXX representative or visit idexx.com/MyLab to locate your local laboratory to learn more.

.....



Choose IDEXX pregnancy tests for proven performance

IDEXX Bovine Pregnancy Test¹

Species	Sample type	Sensitivity	Specificity	Day of detection	Days postpartum
Bovine	Serum	99.3%	93.8%	28	60
	Plasma	99.1%	95.1%	28	60
	Serum (ET)	27/27*	N/A	21 post-embryo transfer	N/A
Water buffalo	Plasma	98%	100%	30	60
Bison	Plasma	90%–94.7% compared to another commercial ELISA		N/A	N/A
Sheep	Serum	100%	100%	35	N/A
Goat	Serum	100%	91.4%	28	N/A

*27 of 27 heifers tested pregnant at 19–21 days post-embryo transfer.

IDEXX Milk Pregnancy Test¹

Species	Sample type	Sensitivity	Specificity	Day of detection	Days postpartum
Bovine	Milk	98.7%	94.4%	28	60
Goat	Milk	94.9%	100%	28	N/A
Sheep	Milk	100%	95%	60	80
Water buffalo	Milk	(4/4)*	100%	29	40

*4 of 4 water buffaloes tested pregnant at 29 days post breeding.

IDEXX Rapid Visual Pregnancy Test¹

Species	Sample type	Sensitivity	Specificity	Day of detection	Days postpartum
Bovine	Serum	98.6%	92.4%	28	60
	Plasma	89.7%	98.9%	28	60
	Whole blood	99%	100%	28	60
Goat	Serum	98.3%	93.7%	28	N/A
Sheep	Serum	99.3%	96.3%	35	N/A
Sheep	Whole blood	100%	98.6%	60	83
Water buffalo	Whole blood	95.9%	94.6%	40	N/A

For more information about these pregnancy tests, refer to the Validation Report for each test.

IDEXX Laboratories, Inc.

One IDEXX Drive
Westbrook, Maine 04092
USA

Tel: +1 207 556 4890 or
+1 800 548 9997
Fax: +1 207 556 4826 or
+1 800 328 5461

IDEXX Europe B.V.

Scorpius 60 Building F
2132 LR Hoofddorp
The Netherlands

Tel: +31 23 558 70 00 or
+800 727 43399
Fax: +31 23 558 72 33

IDEXX Laboratories, Inc.

3F-5 No. 88, Rei Hu Street
Nei Hu District
11494 Taipei
Taiwan

Tel: +886 2 6603 9728
Fax: +886 2 2658 8242

IDEXX Brasil

1478 Av. Brig. Faria Lima
São Paulo, SP
Brasil

Tel: +55 11 3095-5632
Fax: +55 11 3095-5641

Reference:

1. Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA.

© 2017 IDEXX Laboratories, Inc. All rights reserved. • 09-82105-01

All ®/TM marks are owned by IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

The IDEXX Privacy Policy is available at idexx.com.

Test with Confidence™

