Quanti-Tray* Sealer PLUS

User Manual

For technical support, please call:

North/South America: +1 207 556 4496 or +1 800 321 0207

Europe: +00800 4339 9111

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: +81 422 71 5921

Australia: +1300 44 33 99

idexx.com/water

IDEXX

One IDEXX Drive
Westbrook, Maine 04092 USA

Quanti-Tray* Sealer PLUS

User Manual

English

WARNING

RISK OF ELECTRIC SHOCK

Do not remove the outer shell of the sealer.

Do not spill water into the sealer.

When pouring water samples into the Quanti-Tray*, follow instructions in the Quanti-Tray insert.

RISK OF SKIN BURNS

Only trained personnel should clean the inside of the sealer. For cleaning instructions, refer to the "Cleaning and Preventive Maintenance" section of this manual.

Overview

The IDEXX Quanti-Tray* Sealer PLUS is a motor-driven, heated-roller instrument designed to seal water samples inside IDEXX Quanti-Trays. When used with Quanti-Trays and any IDEXX reagent (such as Colilert*, Colilert*-18, Colisure*, Enterolert*, Pseudalert*, or HPC for Quanti-Tray*), the sealer automates the sample handling of bacterial enumeration.

Contents

Each IDEXX Quanti-Tray Sealer PLUS comes with the base sealer unit, regional power cord, input shelf, one Quanti-Tray* Rubber Insert (blue), one Quanti-Tray/2000 Rubber Insert (black), MPN tables, and this user manual. **Use only the blue or black rubber inserts with the Quanti-Tray Sealer PLUS.**



Note: Additional Quanti-Tray Rubber Inserts are available separately.

Technical Specifications of the Instrument

Weight	24 lb	11 kg
Height	13 in.	33 cm
Depth	14 in.	35.5 cm
Width	16 in.	40.6 cm
Operating Temperature	50–90°F	10–32°C
Storage Temperature	–4 to 176°F	–20 to 80°C
Power	100–240 V 50/60 Hz 10 A	

Setting Up the Sealer

- Unpack the contents and save packaging for future shipping.
- Position the sealer on a level surface with the display facing you and with adequate space for trays to be inserted and ejected.
- Ensure that the power switch is turned off.
- Plug the power cord into the sealer and then into a grounded outlet of proper voltage (see product label).
- To connect the **optional** input shelf, center the holes over the corresponding tabs and press down.

Understanding the Sealer Display

Counter: Number of trays that have been run through the sealer.

Amber bars: Sealer is warming up.

Green bar: Sealer has warmed to sealing temperature.

Moon icon ☾: Sleep mode is enabled. See the “Using Sleep Mode” section of this manual.

Alert icon with error code ⚠️: See the “Troubleshooting” section of this manual.

Sealing Quanti-Trays

WARNING: Personal injury, sealer damage, property damage, and/or inaccurate test results may occur if sealer is not used according to the following instructions.

1. Press the power switch to turn on the sealer.
The amber status light illuminates.
2. Allow the sealer to warm up until the amber status light changes to green.
The sealer will not operate until the status light turns green, indicating that the unit has reached sealing temperature.
3. Place a Quanti-Tray or Quanti-Tray/2000 tray filled with sample and IDEXX reagent well-side down onto the corresponding rubber insert, making sure that each well is seated firmly into its corresponding hole.
CAUTION: Do not apply labels of any kind to trays before sealing. Labels could damage the sealer. Use a felt-tip marker, or apply a label after the tray has been sealed.
4. Guide the loaded rubber insert into the input slot of the sealer, with the large cutout facing away from the sealer.
5. Slide the loaded rubber insert into the sealer until the motor grabs the insert and begins to draw it into the sealer.
You can insert multiple rubber inserts consecutively without pausing.
After about 10 seconds, the sealed tray is partially ejected from the rear of the sealer.
Note: If you need to reverse the motor (for example, if a misaligned tray was accidentally fed into the sealer), press and hold the **Reverse** button .
CAUTION: Do not reverse the motor if the rubber insert has been drawn fully into the input slot.
6. Remove the rubber insert and tray from the rear of the sealer.
7. Turn off the sealer or use sleep mode when not in use.

Using Sleep Mode

When sleep mode is enabled, the sealer enters a low power mode after 1 hour of inactivity. When the sealer is reactivated, it warms to sealing temperature faster than from a cold start.

To enable or disable sleep mode: Press the **Moon** button .

When sleep mode is enabled:

- A moon icon  is visible on the display.
- After 1 hour of inactivity, the green status light changes to a pulsing amber light, and the LCD display screen turns off.
- Sleep mode will remain enabled even though the sealer has been turned off and back on. The sealer will “sleep” again after 1 hour of inactivity.

To reactivate the sealer after the LCD display has turned off: Press any button, or insert a Quanti-Tray. The sealer will process the tray once it has warmed to sealing temperature, indicated by a green status light.

To disable sleep mode: Press the Moon button again. The moon icon disappears from the display. Sleep mode can be disabled at any time the sealer is active (i.e., LCD display screen is on).

Troubleshooting

IMPORTANT: If any of the following solutions fails to resolve a problem, contact IDEXX Technical Service at +1 800 321 0207.

Symptom	Cause and/or Recommended Action
 1 Internal Temperature Error	The internal temperature of the sealer is over the allowed temperature. The heater turns off automatically so it can cool to a safe temperature. Overheating can be caused by periods of extremely heavy use or by using the sealer in a high-temperature environment. Wait until the green light appears before continuing.
 2 Temperature Read Error	The roller temperature cannot be accurately read. Power the sealer off and then on.
 3 Access Door Error	The upper access door was opened while the unit was turned on. If the error appears while the door is closed, check that the upper access door and the lower access panel are fully seated.
 4 Voltage Error	The input voltage to the sealer is fluctuating outside the normal range. The problem could be caused by faulty input power. Try another power source. If the problem persists, call IDEXX Technical Service at the number above.
 5 Roller Temperature Error	The temperature of the heated roller is out of the optimal sealing range. The problem could be caused by sealing multiple cold samples one after the other or if the instrument is turned on before it has reached operating temperature. Wait for the green light to appear before continuing to run samples.
 6 Temperature Stability Error	The heated roller temperature jumped between consecutive reads. Check the roller for debris, and clean the roller if needed. See "Cleaning and Preventive Maintenance."
 7 Temperature Error during Seal	The temperature of the heated roller dropped below sealing temperature. After the sealing process is complete, the error disappears, but the power light remains amber until the unit heats to the ready state. The problem could be caused by continuous sealing of cold samples. Wait for the green light to appear before continuing to run samples.
Amber power light does not come on.	<ul style="list-style-type: none"> • Make sure the sealer is plugged into a working outlet with the correct voltage. • Make sure the power switch is on. • Check fuses and replace if necessary. See "Changing Fuses."
Green light takes more than 10 minutes to illuminate.	Call IDEXX Technical Service at +1 800 321 0207 or your local distributor.
Green light goes off during operation.	The temperature of the heated roller is under the optimal sealing range. Wait for the green light to appear before continuing to run samples.
Sealer is on, but motor does not start when a Quanti-Tray is inserted.	<ul style="list-style-type: none"> • Make sure the sealer has reached sealing temperature (green light is on). • Make sure the Quanti-Tray is face down in the rubber insert, with the white side of the tray facing up. • Make sure the rubber insert is pushed firmly into the roller until the motor engages it.
The motor starts when a Quanti-Tray is inserted but does not pull the tray through.	<ul style="list-style-type: none"> • Make sure the Quanti-Tray is correctly seated in the rubber insert and that the insert is pushed firmly into the roller until the motor engages it. • Clean the rubber insert. See the "Cleaning and Preventive Maintenance" section of this manual. • Make sure the insert is completely dry. Wipe it, if necessary.
Quanti-Tray and rubber insert are partially stuck inside the sealer.	Press and hold the Reverse button to reverse the motor until the tray is completely ejected from the input slot. After the rubber insert is removed, make sure that the insert is dry. Wipe it, if necessary. Note: If the tray has gone all the way into the sealer, do not use the Reverse button. Doing so can cause the tray to become lodged inside the sealer. (See the next symptom.)

Symptom	Cause and/or Recommended Action
<p>Quanti-Tray and rubber insert are stuck inside the sealer and are not accessible from the input slot.</p>	<p>Pull the rubber insert out of the exit slot.</p>
<p>Sealer makes a loud hissing noise when sealing, or liquid drips from the bottom of the sealer. Green status light may go out after trays are filled.</p>	<p>The tray was overfilled. The sealer automatically removes excess sample from the tray and discards it into the internal drip tray.</p> <p><i>To prevent the problem:</i></p> <ul style="list-style-type: none"> • Do not fill trays with more than 100 mL of sample. • Avoid excessive water on the surface of the trays. <p><i>To clean the sealer:</i></p> <p>Refer to the “Cleaning and Preventive Maintenance” section of this manual, for instructions on cleaning the inside of the sealer.</p> <p>Important: Do not remove the outer shell of the sealer or tilt the sealer if sample is dripping from the unit. Wait until dripping stops and the sealer is cool. Only trained personnel should clean the sealer.</p>

Cleaning and Preventive Maintenance

Cleaning the Rubber Insert

You can autoclave the rubber insert or clean it with mild detergent, diluted bleach, or isopropyl alcohol, taking the usual precautions when handling such liquids. Be sure to clean the rubber insert if it shows any signs of built-up grime and ensure the rubber insert is dry prior to use.

Cleaning the Outside of the Sealer

Clean the outside of the sealer with a soft, dry cloth or with a soft cloth moistened with water, mild detergent, diluted bleach, or isopropyl alcohol, taking the usual precautions when handling such liquids.

Cleaning the Inside of the Sealer

Routine maintenance and prompt cleaning after a spill will help to maintain the performance of the sealer.

Only qualified personnel should open the access doors to clean the inside of the sealer.

CAUTION: Do not open the access doors of the sealer or tilt the sealer if sample is dripping from the unit. Wait until dripping stops and the sealer is cool.


For troubleshooting or repairs other than cleaning, contact IDEXX Technical Service at +1 800 321 0207.

To clean the inside of the sealer:

1. **Power off the sealer:** Unplug the sealer, and allow it to cool for at least 30 minutes.

2. Open the sealer:

- a) Use a standard flathead screwdriver to unlock the access panel, turning the screw 90° clockwise.
- b) Press in the tabs on both sides of the lower access panel, and then slide the access panel out of the sealer and set aside.

2a.

2b.

- c) Press the round indentation in the upper access door to push the door inward.

2c.

- d) Lift upward to rotate the upper access door fully into the sealer, exposing the rollers. When the upper access door assembly clicks into place, the sealer is ready for cleaning.

2d.

3. Clean the rollers:


CAUTION: Do not touch the upper roller if it is hot.

- a) Soak up excess water from the drip tray and then lift up and remove the drip tray from the bottom of the sealer.

3a.

- b) To remove debris from the rollers, use a soft cloth or sponge soaked with water.
- c) To decontaminate the rollers, use a soft cloth moistened with mild detergent, diluted bleach, or isopropyl alcohol.
- Never use abrasive materials for cleaning.
 - Never use caustic cleaners.
 - Use alcohol only if the sealer is cool.
 - For stubborn deposits, allow time for the cleaner to soak.
 - Do not clean the white insulating layer.

3b.-3c.

- d) Use paper towels or a soft, dry cloth to dry the rollers, interior surfaces, and drip tray.

3d.

4. Close the sealer:


a) Lift the lever slightly to release the upper access door assembly. The upper access door returns to the closed position.

4a.

b) Fully slide the lower access panel back into the sealer until it clips into place.


Note: You may need to lift up on the metal insert platform while sliding the lower access panel into place.

4b.

c) Use a standard flathead screwdriver to lock the panel, turning the screw 90° **counterclockwise**.

Do not overtighten the latch.

4c.

Changing Fuses

Fuses are located in the fuse holder beside the power switch.

If you need to change a fuse, turn off the sealer and unplug the power cord before opening the fuse box.

Use two 10 A, 5 x 20 mm, time-lag fuses. Designed to IEC or local requirements.

LIMITED WARRANTY

IDEXX Laboratories, Inc. ("IDEXX") warrants this product to conform to our published specifications, when stored under appropriate conditions and given normal, proper and intended usage, until the expiration of its stated shelf life, or, if none is stated, for one year from the date of delivery of this product to the original end user purchaser ("Buyer"). IDEXX agrees during the applicable warranty period to replace all nonconforming products within 30 days after date of return to IDEXX and without cost to Buyer. IDEXX shall not have any obligation under this Limited Warranty to make replacements which result, in whole or in part, from catastrophe, fault or negligence of the Buyer, or anyone claiming through or on behalf of the Buyer, or from improper use of the products, or use of the products in a manner for which they were not designed, or by causes external to the products.

Buyer shall notify IDEXX of any products which it believes to be nonconforming during the warranty period. At IDEXX's option, such products shall be returned by Buyer, transportation and insurance prepaid, to IDEXX's designated facility for examination and testing. IDEXX shall repair or replace, within 30 days of receipt by IDEXX, any such product found to be so nonconforming and promptly return such products to Buyer, transportation and insurance prepaid. Should IDEXX's examination and testing not disclose any nonconformity covered by the foregoing warranty, IDEXX shall so advise Buyer and dispose of or return the product in accordance with Buyer's instructions and at Buyer's sole expense.

THE PROVISIONS OF THE FOREGOING LIMITED WARRANTY ARE IN LIEU OF ANY OTHER WARRANTY, WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL (INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE). IDEXX'S LIABILITY ARISING OUT OF THE MANUFACTURE, SALE OR SUPPLYING OF THE PRODUCTS OR THEIR USE OR DISPOSITION, WHETHER BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE, SHALL NOT EXCEED THE ACTUAL PURCHASE PRICE PAID BY BUYER FOR THE PRODUCTS. IN NO EVENT SHALL IDEXX BE LIABLE TO BUYER OR ANY OTHER PERSON OR ENTITY FOR SPECIAL, INCIDENTAL, CONSEQUENTIAL, INDIRECT OR EXEMPLARY DAMAGES (INCLUDING, BUT NOT LIMITED TO, LOSS OF PROFITS OR LOSS OF USE DAMAGES) ARISING OUT OF THE MANUFACTURE, SALE OR SUPPLY OF THE PRODUCTS. THE FOREGOING WARRANTIES EXTEND TO BUYER ONLY AND SHALL NOT BE APPLICABLE TO ANY OTHER PERSON OR ENTITY INCLUDING, WITHOUT LIMITATION, CUSTOMERS OF BUYER.

For technical support, please call:

North/South America: +1 207 556 4496 or +1 800 321 0207

Europe: +00800 4339 9111

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: +81 422 71 5921

Australia: +1300 44 33 99

idexx.com/water

*Quanti-Tray, Colilert, Colisure, Enterolert and Pseudalert are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

© 2015 IDEXX Laboratories, Inc. All rights reserved.

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© 2015 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.