

# IDEXX SNAP tests consistently deliver superior performance



## Introduction

When it comes to infectious disease testing, accurate, timely results are paramount. However, head-to-head comparisons reveal significant differences between the sensitivity and specificity of various in-house tests on the market. These differences can impact your clinical decisions. As this paper demonstrates, IDEXX SNAP® tests consistently outperform other in-house tests in peer-reviewed studies.

## Recent peer-reviewed study results

Antigen/Antibody	IDEXX SNAP® test	Sensitivity	Lateral flow test	Sensitivity	Sample size*	Reference method(s)
Heartworm	SNAP® 4Dx® Plus Test	97.4%	VetScan® FLEX4 Rapid Test <sup>1</sup>	76.9%	232	ELISA <sup>†</sup>
<i>B. burgdorferi</i> (Lyme)	SNAP® 4Dx® Plus Test	95.5%	VetScan® FLEX4 Rapid Test <sup>2</sup>	40.9%	105	Western Blot
<i>Ehrlichia canis</i>	SNAP® 4Dx® Plus Test	97.1%	VetScan® FLEX4 Rapid Test <sup>2</sup>	61.4%	154	ELISA/IFA
<i>Ehrlichia ewingii</i>	SNAP® 4Dx® Plus Test	98.2%	VetScan® FLEX4 Rapid Test <sup>2</sup>	59.3%	163	ELISA
<i>Ehrlichia chaffeensis</i> <sup>‡</sup>	SNAP® 4Dx® Plus Test	64.3%	VetScan® FLEX4 Rapid Test <sup>2</sup>	35.7%	151	ELISA
<i>Anaplasma phagocytophilum</i>	SNAP® 4Dx® Plus Test	84.5%	VetScan® FLEX4 Rapid Test <sup>2</sup>	12.7%	160	IFA
<i>Anaplasma platys</i>	SNAP® 4Dx® Plus Test	83.3%	VetScan® FLEX4 Rapid Test <sup>2</sup>	33.3%	115	ELISA
<i>Giardia</i>	SNAP® <i>Giardia</i> Test	87.1%	VetScan® Canine <i>Giardia</i> Rapid Test <sup>3,4</sup>	70%	101	IFA
<i>Giardia</i>	SNAP® <i>Giardia</i> Test	87.1%	Witness® <i>Giardia</i> Test <sup>4</sup>	73.3%	87	IFA
FeLV	SNAP® FIV/FeLV Combo Test	100%	VetScan® Feline FIV/FeLV Rapid Test <sup>5</sup>	85.6%	146	ELISA
FeLV	SNAP® FIV/FeLV Combo Test	100%	Witness® FeLV-FIV Test <sup>5</sup>	89%	146	ELISA

\*Comparative performance using field samples.

<sup>†</sup>PetChek ELISA with heat treatment.

<sup>‡</sup>The IDEXX SNAP 4Dx Plus Test is not labeled for detection of antibodies to *E. chaffeensis*.

but a peer-reviewed study showed superior sensitivity when compared to the VetScan FLEX4 Rapid Test.

## SNAP tests leverage reference laboratory quality

IDEXX scientists created the SNAP® assay by combining advanced technology used in reference laboratories with well-defined diagnostic markers. While the platform of SNAP tests continues to expand, each test delivers on the promise of reference-laboratory quality with pet-side results.

## Three unique steps enhance the sensitivity and specificity of SNAP tests to help you make the right decision for your patient

### Bidirectional flow enhances sensitivity


The patient's sample is mixed with a liquid specifically formulated for each test. When added to the SNAP® device, this mixture flows over the capture reagents on the test spots toward the activation circle. Activating the device draws the sample mixture back across the capture reagents a second time. This provides an additional opportunity for binding on the test spots, which increases the sensitivity of results.

### Integrated wash step maximizes specificity

Activating the SNAP device also triggers a wash buffer that removes unbound sample (blood cells, fecal material) from the surface of the device. This enhances the specificity of the test and produces a clean white background that makes it easy to interpret results.

### Amplification reveals low-level positives

Once the sample mixture has bound to the test spots, it undergoes an amplification process to maximize the signal, further increasing the sensitivity of the test.



## Summary and conclusions

The unique SNAP® platform, with its patented bidirectional flow, integrated wash step, and amplification, delivers the superior sensitivity and specificity of reference laboratory-quality testing to the pet-side assay. For more than 25 years, veterinarians, scientists, and researchers have relied on SNAP® tests to establish the prevalence of infectious diseases worldwide. Their work has validated the accuracy and consistent performance of SNAP tests in over 100 peer-reviewed studies.<sup>6</sup> Today, more than 20 million SNAP tests are run each year, making it one of the most trusted and commonly used in-clinic diagnostic platforms in the world.<sup>6</sup>

For more information, visit [idexx.com/snap](https://idexx.com/snap)

## References

1. Burton KW, Liu J, Drexel JP, Pulaski C, Heartsill S, Dupuy C. Comparative evaluation of field samples using 2 in-clinic assays for heartworm antigen detection in dogs [published online ahead of print June 17, 2020]. *Vet Parasitol.* 2020;283:109163. doi:10.1016/j.vetpar.2020.109163
2. Liu J, Drexel J, Andrews B, Eberts M, Breitschwerdt E, Chandrashekar R. Comparative Evaluation of 2 In-Clinic Assays for Vector-Borne Disease Testing in Dogs. *Top Companion Anim Med.* 2018;33(4):114-118. doi:10.1053/j.tcam.2018.09.003
3. Bowman DD, Liotta JL, Mizhquiri-Barbecho JS, Rishniw M, Simpson K. Performance of four in-clinic tests for *Giardia* antigen compared with direct immunofluorescence assay. *J Vet Intern Med.* 2017;31(5):1590. doi:10.1111/jvim.14784
4. Barbecho JM, Bowman DD, Liotta JL. Comparative performance of reference laboratory tests and in-clinic tests for *Giardia* in canine feces. *Parasit Vectors.* 2018;11(1):444. doi:10.1186/s13071-018-2990-6
5. Levy JK, Crawford PC, Tucker SJ. Performance of 4 point-of-care screening tests for feline leukemia virus and feline immunodeficiency virus. *J Vet Intern Med.* 2017;31(2):521-526. doi:10.1111/jvim.14648
6. Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA.