English version

SNAP® Foal IgG
The SNAP® Foal IgG Test is an enzyme immunoassay for semi-quantitative detection of immunoglobulin G (IgG) in equine serum, plasma or whole blood.

Storage
SNAP devices and test reagents are stable until the expiration dates when stored at 2-8°C. Do not allow SNAP devices or test reagents to freeze.

Precautions and Warnings
Use a separate sample diluent bottle and sample loop for each test. For whole-blood samples, use two separate sample loops for each test. Do not reuse sample loops.

Instructions for Use
Kit reagents and devices must be at 18-25°C when used. Whole blood, serum or plasma can be used. Whole blood must be anticoagulated with EDTA or heparin, and can be used either fresh or refrigerated at 2-8°C for up to one week.

Test Procedure
IMPORTANT: Do not depress the activator until indicated. Remove the cap and dropper tip from the sample diluent bottle. Gently mix whole-blood samples by inverting.

For whole-blood samples, we recommend immersing the loop tip in the sample that remains in the cap of the sample collection container. Immerse the loop tip only. Do not immerse the loop handle in the sample.

Transfer the filled loop by immersing and twirling the loop tip in the bottle of sample diluent.

Firmly seat the dropper tip on the sample diluent bottle. Mix thoroughly by inverting five times. Hold the sample diluent bottle vertically and dispose of the first 5-10 drops from the bottle.

Place the SNAP device on a flat surface. With the bottle tip 1/2 to 1 inch directly above the SNAP device, carefully apply one drop of diluted sample directly onto the sample spot in the result window.

Remove the cap from the conjugate bottle and pour its contents into the sample well of the SNAP device. (Some of the contents will remain in the conjugate bottle.) The sample will flow across the result window, reaching the activation circle in 30 to 90 seconds.

Watch the device carefully for color in the activation circle. When color FIRST appears in the activation circle, push the activator firmly until it is flush with the device body.

Wait 7 minutes. Visually read the test result.

Interpreting the Test Result
To determine the test result, read the reaction spots in the result window. To determine the level of IgG in the test sample, compare the color of the sample spot to the color of the two calibrator spots.

Color calibration chart showing reaction spots for IgG levels: <400 mg/dL, Less than 400 mg/dL IgG, Approximately 400 mg/dL IgG, Greater than 800 mg/dL IgG, 400 to 800 mg/dL IgG.


Invalid Test Results
Improper Sample Application
If the drop of diluted sample applied to the SNAP device does not wet the sample spot completely, nonuniform color development on the sample spot can occur.

Background
If the SNAP device is activated too late (the sample is allowed to flow past the activation circle), high background color may result. If colored background obscures the test result, repeat the test.

No Color Development
If no color develops on the calibrator spots, repeat the test. If no color develops on the calibrator spots for a second test device, contact IDEXX Technical Support.

Color Does Not Appear in the Activation Circle
If color does not appear in the activation circle within 90 seconds after the conjugate is added to the SNAP device, repeat the test. If this also happens for a second test device, contact IDEXX Technical Support.

Symbol descriptions
Use by date, Lot, Serial number, Temperature limitation, Manufacturer, Catalog number.

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*SNAP is a trademark or a registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.