Canine Leptospira Antibody Test Kit

For veterinary use only.

SNAP* Lepto

In vitro diagnostic for the detection of anti-Leptospira antibodies to the serovars Grippotyphosa, Canicola, Pomona, and Icterohaemorrhagiae in canine serum.

Precautions and Warnings

- All wastes should be properly decontaminated prior to disposal.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- This test kit cannot differentiate between antibodies to pathogenic and nonpathogenic serovars (Autumnalis and Bratislava).
- This test kit cannot differentiate between antibodies produced as a result of a natural infection with the organism and those produced by vaccination.
- WARNING:** Conjugate - H316/P332+P313/EUH208. Causes mild skin irritation. If skin irritation occurs: Get medical advice/attention. Contains Kathon. May produce an allergic reaction.

Storage

- Store at 2–8°C until expiration date.
- All components must be at room temperature (18–25°C) before running the test.

Kit Components


Item	Reagents	Quantity
1	1 bottle of <i>Leptospira</i> Antigen:HRPO Conjugate (Preserved with gentamicin and Kathon)	3.8 mL
2	SNAP Device	5 or 10
Reagents contained in each device:		
	Wash Solution (preserved with Kathon)	0.4 mL
	Substrate Solution	0.6 mL
Other Components: transfer pipettes, sample tubes and reagent rack		

Sample Information

- Samples must be at room temperature 18–25°C before beginning the test procedure.
- Serum either fresh or stored at 2–8°C for up to one week, can be used.
- For longer storage, serum can be frozen (-20°C or colder) and then reconstituted before use.
- Hemolyzed or lipemic samples will not affect test results.

Test Procedure

- Allow all components to equilibrate at room temperature 18–25°C for 30 minutes before use. **Do not heat.**

Note: Some samples may not flow to the activation circle within 60 seconds, and therefore, the circle may not turn color. In this case, press the activator after the sample has flowed across the result window.


- Read test result at **10 minutes**.

Note: The positive control may develop sooner, but results are not complete until 10 minutes.

Interpreting the Test Results

Positive Result

Any color development in the sample spots indicates the presence of anti-*Leptospira* antibody in the sample.

Note: A positive result indicates presence of antibodies to any of the following serovars: Grippotyphosa, Canicola, Pomona, and Icterohaemorrhagiae. The sample spot cannot differentiate between the four serovars.

Negative Result


Only positive control spot develops color.

Invalid Results

- Background** — If the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No Color Development** — If the positive control does not develop color, repeat the test.

Performance Summary by MAT Titer

Symbol Descriptions

Use by date	ECREP	Authorized Representative in the European Community
LOT Batch Code (Lot)		
SN Serial number	IVD	Consult instructions for use
Temperature limitation	IVD	In vitro diagnostic
Manufacturer	REF	Date of manufacture
Catalog Number		

IDEXX Technical Support

USA/Canada 1 800 248 2483 • Europe 00800 1234 3399
Australia 1300 44 33 99
idexx.com

U.S. Vet. License No. 313
Product Code: 5052.50

*SNAP is a trademark or registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2015 IDEXX Laboratories, Inc. All rights reserved. 06-00000776-03

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