## Using the **StatSpin**\* Centrifuge Shield

Please retain this notification with your operator's manual. For an updated version, visit [idexx.com/statspin](http://idexx.com/statspin).

### About StatSpin Centrifuge Alarms

The system beeps continuously if one of the following conditions is present:

- The centrifuge temperature exceeds the operating temperature.†
- There is a short circuit of the motor drives, fan, or solenoid.
- There is a reduction in the availability of electrical power.
- There have been short-term power failures.

To stop the beeping, press the **Stop** button.

†If the centrifuge temperature exceeds the operating temperature, you must wait 15 minutes before pressing the **Stop** button.

### Intended Use

The shield is for use with the following StatSpin\* Centrifuge models:

- M900-10
- M900-22
- V900-22
- M900-12
- M901-22
- V901-22

### Installing the Shield


A shield should always be used with the RT12 rotor and should be replaced when the rotor is replaced.

1. Ensure the RT12 rotor is properly seated.
2. Press the shield into the bowl until it makes full contact with the bottom of the bowl.

**Note:** To remove the shield, pull up on the center edge of the shield.

### Shield Care and Maintenance

- Do not operate the centrifuge without a shield in place when using an RT12 rotor.
- Do not use a shield with the RH12 microhematocrit rotor.
- Always install a new shield when installing a new RT12 rotor.
- Always replace the shield when replacing the RT12 rotor or when any crack or visible wear occurs. Refer to the operator's manual for rotor maintenance frequency.
- To clean, remove the shield from the centrifuge and clean it with a damp cloth and mild detergent.

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