

**VMX**  
2025  
VETERINARY MEETING & EXPO

# Festival of the HeARTS



Presented By:

**NAVCH**  
YOUR VETERINARY COMMUNITY

JANUARY 25-29

ORLANDO, FLORIDA

NAVCOM.COM

**“Test something, no, test  
EVERYTHING....STAT!”**

**Choosing and interpreting useful  
diagnostics in emergency patients**

Nancy A. Sanders, DVM  
DACVIM (SAIM), DACVECC

## **Conflict of Interest Disclosure:**

Nancy Sanders is a full-time IDEXX employee

The information contained herein is intended to provide general guidance only. Diagnosis, treatment, and monitoring should be patient specific and is the responsibility of the veterinarian providing primary care.

# Demonstrate prioritization of ER diagnostics by use of:

- 1 | Assessment of patient status
- 2 | Your raw talents
- 3 | Determining patient and pet guardian needs
- 4 | Your experiences
- 5 | Case examples

Before we talk about which  
tests...

Let's talk about when



# Never delay necessary treatment to run diagnostics

- + Stabilize respiratory status
  - + Ensure patent airway
  - + Oxygen by least stressful method
  - + Thoracocentesis
  - + Sedation if indicated
- + Support cardiovascular system
  - + Intravenous or intraosseous fluids
  - + Medications
  - + Stop bleeding, blood products
  - + Pericardiocentesis
- + Cover pain, anxiety, stress
  - + Narcotic analgesics
  - + Acepromazine (low dose, cautious)

**A**irway

**B**reathing

**C**ardiovascular

# Diagnostics for who?

Patient

Pet guardian

Your ER's protocols

The referring DVM/practice

Factors that  
influence diagnostic  
prioritization

# Define the purpose of your emergency room


## For example...

- + Stabilize and send home or to referring vet
  - + “Treat and street”
- + Do everything
  - + Full diagnosis
  - + Primary vet
- + Referred for specific diagnostic or overnight monitoring “only”
- + Urgent vs. emergent care

# The balancing act



# Choosing diagnostics for the patient



What's the client complaint?  
What's the patient's main problem?  
What's the emergency?

FOCUS!!!



# What constitutes an emergency (a.k.a. what kills rapidly?)

## The BIG 3

### Cardiovascular

- CHF
- Arrhythmias
- Thromboembolism
- Hemorrhage
- Shock (any form)
- Hypotension
- Hypertension

### Neurologic

- Seizures
- Trauma
- Stroke
- Increased intracranial pressure
- Hypertension
- Hypotension

### Respiratory

- Hypoxia
  - Pulmonary infiltrates
  - Pleural space disorders
  - Thromboembolism
  - Severe hypoventilation
- Hypercapnia
  - Hypoventilation of any cause



# Which diagnostics and when? (Examples only!)

## Immediate – 1<sup>st</sup> tier

- Quick physical/ patient assessment
- PCV/TP
- Blood glucose (BG)
- Serum electrolytes
- ECG
- Systolic blood pressure (SBP)
- Pulse oximetry
- Point of care ultrasound (POCUS)


## ASAP – 2<sup>nd</sup> tier


- CBC/chemistry/urinalysis
- Lactate
- Arterial or venous blood gas
- Coagulation profile
- SNAP® proBNP
- Imaging
- Fluid samples/analysis/cytology

## When appropriate - 3<sup>rd</sup> tier


- Quant PL/ Snap® cPL/fPL
- Infectious disease screening
- C-reactive protein
- Cystatin B
- Blood type/cross-match
- Blood smear
- Slide agglutination
- Serum cortisol
- Advanced imaging
- ProBNP

**\*\*Extra samples\*\***







- Primary central: CSF, advanced imaging
- Secondary central: BP, coag panel
- Metabolic: BG,  $iCa^{++}$ , chemistry, bile acids, UA
- Toxic: CBC/chemistry/UA



- Upper airway: Oral exam, imaging, venous blood gas (VBG)
- Lower airway: Chest rads, ABG, lung sampling
- Pleural space: TFAST, pleural space sampling/evacuation
- Really pale patient: PCV/TP, SBP
- Suspect cardiac: Rads, TFAST, ECG, proBNP



- Vomiting: Abdominal imaging, labs
- Regurgitation: Chest rads
- Ptyalism: Oral exam vs. GI work-up




- Large bowel diarrhea: rectal, fecal, abdominal imaging
- Small bowel diarrhea: blood work including GI panel, fecal, abdominal imaging

# When considering diagnostics, a little forward-thinking can go a LONG way

- What question(s) am I trying to answer?
- What do results of this test tell me?
- Will results of this test change what I do?
- What might happen to the patient if treatment is withheld?
- What might happen if treatment is administered without full investigation?
- How might current treatments interfere with future diagnostics?
- What future/additional samples might I need?
- If saving samples, how should they be stored?

# General “minimum database”




# Hematology consists of 3 parts


## CBC

Hematology		5/13/20	11:25 AM
<a href="#">Click to view Differentials</a>			
■ RBC	6.33	6.54 - 12.20 Mj/L	
■ Hematocrit	29.2	30.3 - 52.3 %	
■ Hemoglobin	9.7	9.8 - 16.2 g/dL	
■ MCV	46.1	35.9 - 53.1 fL	
■ MCH	15.3	11.8 - 17.3 pg	
■ MCHC	33.2	28.1 - 35.8 g/dL	
■ RDW	20.9	15.0 - 27.0 %	
■ % Reticulocyte	0.1	%	
■ Reticulocytes	5.7	3.0 - 50.0 Kj/L	
■ Reticulocyte Hemoglobin	15.2	13.2 - 20.8 pg	
■ WBC	8.91	2.87 - 17.02 Kj/L	
■ % Neutrophils	56.9	%	
■ % Lymphocytes	33.4	%	
■ % Monocytes	5.6	%	
■ % Eosinophils	3.4	%	
■ % Basophils	0.7	%	
■ Neutrophils	5.07	2.30 - 10.29 Kj/L	
■ Lymphocytes	2.98	0.92 - 6.88 Kj/L	
■ Monocytes	0.50	0.05 - 0.67 Kj/L	
■ Eosinophils	0.30	0.17 - 1.57 Kj/L	
■ Basophils	0.06	0.01 - 0.26 Kj/L	
■ Platelets	368	151 - 600 Kj/L	
■ MPV	16.7	11.4 - 21.6 fL	
■ Plateletcrit	0.61	0.17 - 0.86 %	

## Dot plots/graphics



## Blood smear





# The three major sections of the CBC

## Erythrogram



### Hematocrit, RBC

Anemia  
Erythrocytosis


### Anemia Classification

Reticulocytes  
Regenerative  
Non-regenerative

### RBC Indices

MCV  
RDW  
Retic-Hemoglobin

## Leukogram



### Total White Blood Count

Cell Distribution

### Cell Descriptions

Left Shift  
Toxic Changes  
Precursor

## Thrombogram



### Total Platelet Count

### Plateletcrit

### Platelet indices

MPV

### Morphology

Megakaryocytes  
Clumping

# How do I analyze the chemistry panel?

Chemistry		Yesterday 12:56 PM		
Glucose	73	72 - 175 mg/dL	<div><div></div><div></div><div></div></div>	
IDEXX SDMA	a. 8	0 - 14 µg/dL	<div><div></div><div></div><div></div></div>	
Creatinine	1.2	0.9 - 2.3 mg/dL	<div><div></div><div></div><div></div></div>	
BUN	35	16 - 37 mg/dL	<div><div></div><div></div><div></div></div>	
BUN: Creatinine Ratio	29.2			
Phosphorus	4.7	2.9 - 6.3 mg/dL	<div><div></div><div></div><div></div></div>	
Calcium	9.6	8.2 - 11.2 mg/dL	<div><div></div><div></div><div></div></div>	
Sodium	150	147 - 157 mmol/L	<div><div></div><div></div><div></div></div>	
Potassium	4.3	3.7 - 5.2 mmol/L	<div><div></div><div></div><div></div></div>	
Na: K Ratio	35	29 - 42		
Chloride	116	114 - 126 mmol/L	<div><div></div><div></div><div></div></div>	
TCO2 (Bicarbonate)	20	12 - 22 mmol/L	<div><div></div><div></div><div></div></div>	
Anion Gap	18	12 - 25 mmol/L	<div><div></div><div></div><div></div></div>	
Total Protein	7.9	6.3 - 8.8 g/dL	<div><div></div><div></div><div></div></div>	
Albumin	3.4	2.6 - 3.9 g/dL	<div><div></div><div></div><div></div></div>	
Globulin	4.5	3.0 - 5.9 g/dL	<div><div></div><div></div><div></div></div>	
Albumin: Globulin Ratio	0.8	0.5 - 1.2	<div><div></div><div></div><div></div></div>	
ALT	45	27 - 158 U/L	<div><div></div><div></div><div></div></div>	
AST	28	16 - 67 U/L	<div><div></div><div></div><div></div></div>	
ALP	29	12 - 59 U/L	<div><div></div><div></div><div></div></div>	
GGT	<1	0 - 6 U/L	<div><div></div><div></div><div></div></div>	
Bilirubin - Total	0.1	0.0 - 0.3 mg/dL	<div><div></div><div></div><div></div></div>	
Bilirubin - Unconjugated	0.0	0.0 - 0.2 mg/dL	<div><div></div><div></div><div></div></div>	
Bilirubin - Conjugated	<0.1	0.0 - 0.2 mg/dL	<div><div></div><div></div><div></div></div>	
Cholesterol	242	91 - 305 mg/dL	<div><div></div><div></div><div></div></div>	
Amylase	981	623 - 2,239 U/L	<div><div></div><div></div><div></div></div>	
Lipase	b. 11	0 - 45 U/L	<div><div></div><div></div><div></div></div>	
Creatine Kinase	223	64 - 440 U/L	<div><div></div><div></div><div></div></div>	
Hemolysis Index	c. N			
Lipemia Index	d. 1+			

Metabolism, endocrine

Kidney

Mineral/Bone

Electrolytes

Acid-base

Proteins

Hepatobiliary

Pancreas, GI

Muscle

Sample Quality

Grouped per basic categories

# Kidney assessment is complex and is a combination of categories



## Kidney Disease

	Metabolism
	Kidney
	Mineral/Bone
	Electrolytes
	Acid-base
	Proteins
	Hepatobiliary
	Pancreas, GI
	Muscle
	Sample Quality

Glucose
Kidney
Mineral
Electrolytes
Acid-base
Proteins
Cholesterol

Chemistry				Yesterday 12:56 PM
Glucose	73	72 - 175 mg/dL		
IDEXX SDMA	a. 8	0 - 14 µg/dL		
Creatinine	1.2	0.9 - 2.3 mg/dL		
BUN	35	16 - 37 mg/dL		
BUN: Creatinine Ratio	29.2			
Phosphorus	4.7	2.9 - 6.3 mg/dL		
Calcium	9.6	8.2 - 11.2 mg/dL		
Sodium	150	147 - 157 mmol/L		
Potassium	4.3	3.7 - 5.2 mmol/L		
Na: K Ratio	35	29 - 42		
Chloride	116	114 - 126 mmol/L		
TCO2 (Bicarbonate)	20	12 - 22 mmol/L		
Anion Gap	18	12 - 25 mmol/L		
Total Protein	7.9	6.3 - 8.8 g/dL		
Albumin	3.4	2.6 - 3.9 g/dL		
Globulin	4.5	3.0 - 5.9 g/dL		
Albumin: Globulin Ratio	0.8	0.5 - 1.2		
ALT	45	27 - 158 U/L		
AST	28	16 - 67 U/L		
ALP	29	12 - 59 U/L		
GGT	<1	0 - 6 U/L		
Bilirubin - Total	0.1	0.0 - 0.3 mg/dL		
Bilirubin - Unconjugated	0.0	0.0 - 0.2 mg/dL		
Bilirubin - Conjugated	<0.1	0.0 - 0.2 mg/dL		
Cholesterol	242	91 - 305 mg/dL		
Amylase	981	623 - 2,239 U/L		
Lipase	b. 11	0 - 45 U/L		
Creatine Kinase	223	64 - 440 U/L		
Hemolysis Index	c. N			
Lipemia Index	d. 1+			

# Liver assessment is complex and a combination of categories

## Chemistry Yesterday 12:56 PM

Glucose	73	72 - 175 mg/dL	
IDEXX SDMA	a. 8	0 - 14 µg/dL	
Creatinine	1.2	0.9 - 2.3 mg/dL	
BUN	35	16 - 37 mg/dL	
BUN: Creatinine Ratio	29.2		
Phosphorus	4.7	2.9 - 6.3 mg/dL	
Calcium	9.6	8.2 - 11.2 mg/dL	
Sodium	150	147 - 157 mmol/L	
Potassium	4.3	3.7 - 5.2 mmol/L	
Na: K Ratio	35	29 - 42	
Chloride	116	114 - 126 mmol/L	
TCO2 (Bicarbonate)	20	12 - 22 mmol/L	
Anion Gap	18	12 - 25 mmol/L	
Total Protein	7.9	6.3 - 8.8 g/dL	
Albumin	3.4	2.6 - 3.9 g/dL	
Globulin	4.5	3.0 - 5.9 g/dL	
Albumin: Globulin Ratio	0.8	0.5 - 1.2	
ALT	45	27 - 158 U/L	
AST	28	16 - 67 U/L	
ALP	29	12 - 59 U/L	
GGT	<1	0 - 6 U/L	
Bilirubin - Total	0.1	0.0 - 0.3 mg/dL	
Bilirubin - Unconjugated	0.0	0.0 - 0.2 mg/dL	
Bilirubin - Conjugated	<0.1	0.0 - 0.2 mg/dL	
Cholesterol	242	91 - 305 mg/dL	
Amylase	981	623 - 2,239 U/L	
Lipase	b. 11	0 - 45 U/L	
Creatine Kinase	223	64 - 440 U/L	
Hemolysis Index	c. N		
Lipemia Index	d. 1+		


Metabolism
Kidney
Mineral/Bone
Electrolytes
Acid-base
Proteins
Hepatobiliary
Pancreas, GI
Muscle
Sample Quality



## Liver Disease

Glucose
Kidney
Electrolytes
Proteins
Hepatobiliary
Pancreas, GI
Sample quality

# GI assessment is complex and often a combination of categories



## GI & Pancreatic Disease

Glucose

Kidney

Mineral/Bone

Electrolytes

Acid-base

Proteins

Hepatobiliary

Pancreas, GI

Sample quality

IDEXX

Metabolism

Kidney

Mineral/Bone

Electrolytes

Acid-base

Proteins

Hepatobiliary

Pancreas, GI

Muscle

Sample Quality

### Chemistry

Yesterday  
12:56 PM

Glucose	73	72 - 175 mg/dL	<div><div></div></div>
IDEXX SDMA	a. 8	0 - 14 µg/dL	<div><div></div></div>
Creatinine	1.2	0.9 - 2.3 mg/dL	<div><div></div></div>
BUN	35	16 - 37 mg/dL	<div><div></div></div>
BUN: Creatinine Ratio	29.2		
Phosphorus	4.7	2.9 - 6.3 mg/dL	<div><div></div></div>
Calcium	9.6	8.2 - 11.2 mg/dL	<div><div></div></div>
Sodium	150	147 - 157 mmol/L	<div><div></div></div>
Potassium	4.3	3.7 - 5.2 mmol/L	<div><div></div></div>
Na: K Ratio	35	29 - 42	<div><div></div></div>
Chloride	116	114 - 126 mmol/L	<div><div></div></div>
TCO2 (Bicarbonate)	20	12 - 22 mmol/L	<div><div></div></div>
Anion Gap	18	12 - 25 mmol/L	<div><div></div></div>
Total Protein	7.9	6.3 - 8.8 g/dL	<div><div></div></div>
Albumin	3.4	2.6 - 3.9 g/dL	<div><div></div></div>
Globulin	4.5	3.0 - 5.9 g/dL	<div><div></div></div>
Albumin: Globulin Ratio	0.8	0.5 - 1.2	<div><div></div></div>
ALT	45	27 - 158 U/L	<div><div></div></div>
AST	28	16 - 67 U/L	<div><div></div></div>
ALP	29	12 - 59 U/L	<div><div></div></div>
GGT	<1	0 - 6 U/L	<div><div></div></div>
Bilirubin - Total	0.1	0.0 - 0.3 mg/dL	<div><div></div></div>
Bilirubin - Unconjugated	0.0	0.0 - 0.2 mg/dL	<div><div></div></div>
Bilirubin - Conjugated	<0.1	0.0 - 0.2 mg/dL	<div><div></div></div>
Cholesterol	242	91 - 305 mg/dL	<div><div></div></div>
Amylase	981	623 - 2,239 U/L	<div><div></div></div>
Lipase	b. 11	0 - 45 U/L	<div><div></div></div>
Creatine Kinase	223	64 - 440 U/L	<div><div></div></div>
Hemolysis Index	c. N		
Lipemia Index	d. 1+		

# Endocrine assessment is complex and often a combination of categories



## Endocrine Disease

Glucose

Kidney

Mineral/Bone

Electrolytes

Acid-base

Proteins

Hepatobiliary

Muscle

Sample quality

IDEXX

Metabolism

Kidney

Mineral/Bone

Electrolytes

Acid-base

Proteins

Hepatobiliary

Pancreas, GI

Muscle

Sample Quality

Chemistry		Yesterday 12:56 PM	
Glucose	73	72 - 175 mg/dL	
IDEXX SDMA	a. 8	0 - 14 µg/dL	
Creatinine	1.2	0.9 - 2.3 mg/dL	
BUN	35	16 - 37 mg/dL	
BUN: Creatinine Ratio	29.2		
Phosphorus	4.7	2.9 - 6.3 mg/dL	
Calcium	9.6	8.2 - 11.2 mg/dL	
Sodium	150	147 - 157 mmol/L	
Potassium	4.3	3.7 - 5.2 mmol/L	
Na: K Ratio	35	29 - 42	
Chloride	116	114 - 126 mmol/L	
TCO2 (Bicarbonate)	20	12 - 22 mmol/L	
Anion Gap	18	12 - 25 mmol/L	
Total Protein	7.9	6.3 - 8.8 g/dL	
Albumin	3.4	2.6 - 3.9 g/dL	
Globulin	4.5	3.0 - 5.9 g/dL	
Albumin: Globulin Ratio	0.8	0.5 - 1.2	
ALT	45	27 - 158 U/L	
AST	28	16 - 67 U/L	
ALP	29	12 - 59 U/L	
GGT	<1	0 - 6 U/L	
Bilirubin - Total	0.1	0.0 - 0.3 mg/dL	
Bilirubin - Unconjugated	0.0	0.0 - 0.2 mg/dL	
Bilirubin - Conjugated	<0.1	0.0 - 0.2 mg/dL	
Cholesterol	242	91 - 305 mg/dL	
Amylase	981	623 - 2,239 U/L	
Lipase	b. 11	0 - 45 U/L	
Creatine Kinase	223	64 - 440 U/L	
Hemolysis Index	c. N		
Lipemia Index	d. 1+		

A close-up, low-angle shot of a diverse group of young people at what appears to be a concert or a large gathering. In the foreground, a young woman with dark skin and curly hair is shown in profile, her mouth wide open in a shout or cheer, with her eyes closed. To her left, another person's face is partially visible, also smiling and cheering. In the background, a young man with blonde hair is seen from the side, his head tilted back and arms raised in the air. The lighting is warm and golden, suggesting an outdoor setting during sunset or sunrise. The overall mood is one of joy, excitement, and collective energy.

Cases!!!

“We’re  
not  
done”



My pet collapsed.



# Jasper

- 15-year-old MC mixed-breed dog
- Collapsed this am
- Similar episode a few weeks ago, spontaneously resolved




# Triage

## + Quick PE

- + TPR = 99 F, 160 bpm, 60 bpm
- + Sternal recumbency, QAR, mentally appropriate
- + Pale, moist MM, bounding pulses
- + 2/6 systolic heart murmur, normal lung auscultation
- + Questionable abdominal fluid wave

## + Quick Assessment

- + Tachycardic, tachypneic
- + Shocky, unstable
- + Ddx include: anemia, hemorrhage/hemoabdomen



# What diagnostics?

- + FIRST STABILIZE!!!

- + Fluids/volume replacement
- + Blood transfusion

- + 1<sup>st</sup> tier diagnostics: PCV/TP, SBP, AFAST vs. abdominal rads (vs. blind abdominocentesis)

- + 2<sup>nd</sup> tier diagnostics:

- + 3<sup>rd</sup> tier diagnostics:

## Regenerative anemia...

(Note normal RBC indices)

Test	Results	Reference Interval	LOW	NORMAL	HIGH
ProCyte Dx (May 13, 2019 2:41 AM)					
RBC	2.28 M/ $\mu$ L	5.65 - 8.87	LOW		
HCT	15.9 %	37.3 - 61.7	LOW		
HGB	5.1 g/dL	13.1 - 20.5	LOW		
MCV	69.7 fL	61.6 - 73.5			
MCH	22.4 pg	21.2 - 25.9			
MCHC	32.1 g/dL	32.0 - 37.9			
RDW	18.7 %	13.6 - 21.7			
%RETIC	18.2 %				
RETIC	153.7 K/ $\mu$ L	10.0 - 110.0	HIGH		


PCV 19%/TS 5.2

# Focus

- + What's the client complaint? Collapsing
- + What's the main problem? Profound anemia, regenerative; hypoproteinemia
- + What's the emergency? Anemia, hypovolemia, hypovolemic shock, possible surgical emergency
- + Do these all fit together; are we addressing all concerns?
  
- + What now?
  - + Treat shock and anemia: transfusion and fluid therapy

# What now (continued)?

- + 2<sup>nd</sup> tier diagnostics: determine cause of regenerative anemia > hemorrhage or hemolysis; assess any consequences of shock (CBC, chemistry, UA, lactate, cystatin B)
- + 3<sup>rd</sup> tier diagnostics: If hemorrhage, determine surgical vs. non-surgical (imaging), coags
- + Assess response to ER treatment
- + IF surgical case, determine if surgery an option: serial SBP/PCV/TP/lactate, ECG, thoracic radiographs/met check



# Next diagnostic steps

- + Additional imaging vs. blind abdominocentesis
- + Fluid analysis: PCV/TP, see if clots, +/- cytology or other diagnostics

Abd. fluid  
PCV 29%/TS 4.7  
Does not clot



Not all causes of hemorrhage are surgical...

IDEXX VetConnect<sup>+</sup> PLUS

Home | Lab Services | Telemedicine | Reference

ABIGAIL

2012 OCT 11 OCT 11 OCT 8 OCT 7 OCT 6 OCT 5 OCT 5 OCT 5 OCT 5



Trending


10/5/12 @ 3:11 pm

Coag Dx Analyzer  
[Show Details](#)

Citrated Prothrombin Time (PT)	> 100.0	11.0 - 17.0 seconds	<input type="text"/>
Citrated Partial Thromboplastin Time (PTT)	> 300.0	72.0 - 102.0 seconds	<input type="text"/>

Don't forget your coags!









# What IF...

Test	Results	Reference Interval	LOW	NORMAL	HIGH
ProCyte Dx (May 13, 2019 2:41 AM)					
RBC	2.28 M/ $\mu$ L	5.65 - 8.87	LOW		
HCT	15.9 %	37.3 - 61.7	LOW		
HGB	5.1 g/dL	13.1 - 20.5	LOW		
MCV	69.7 fL	61.6 - 73.5			
MCH	22.4 pg	21.2 - 25.9			
MCHC	32.1 g/dL	32.0 - 37.9			
RDW	18.7 %	13.6 - 21.7			
%RETIC	18.2 %				
RETIC	153.7 K/ $\mu$ L	10.0 - 110.0	HIGH		

PCV 19%/TS 8.2

# Reticulocytosis:

## From where can we get more information?



### Dot Plot Legend


- |            |              |
|------------|--------------|
| ● RBC_FRAG | ● QUALIBEADS |
| ● RBC      | ● DOUBLET    |
| ● RETICS   | ● PLT        |

# If hemolysis is suspected, diagnostics should include:

- + Blood smear review
- + Slide agglutination
- + Chest and abdominal imaging
- + +/- Infectious disease testing
- + +/- Coomb's test

Rouleaux

Agglutination



# Options for rapid evaluation of blood smears are numerous

- + In-house microscopy, self-reviewed
- + Digital cytology, remote evaluation
- + In-house automated analyzer (artificial intelligence)



# Options for screening for infectious diseases can be overwhelming

## Antibody (Ab), serology



- **+** Indicates:
  - Exposure
- **Limitations**
  - Current infx not confirmed
  - May need paired titer

## Antigen (Ag)




- **+** Indicates:
  - Organism present
- **Limitations**
  - Does not confirm organism viability

## PCR



- **+** Indicates:
  - Organism present
- **Limitations**
  - Does not confirm organism viability
  - TAT

## Cytology Histopathology




- **+** Indicates:
  - Organism present
- **Limitations**
  - Specific ID
  - \$

## Culture



- **+** Indicates
  - **Viable** organism
- **Limitations**
  - Sample handling
  - \$
  - TAT


## Other advanced diagnostics





- **+** Indicates:
  - Organism present
- **Limitations**
  - TAT
  - \$
  - Availability

More than one test method and/or repeat testing **maybe be** necessary for a diagnosis

Clinical  
picture  
influences  
test choice



Not all hemolytic anemias are immune-mediated...



# Hemolysis

## Hemolytic anemia (HA) (NOT immune-mediated)

- + Intravascular (hemoglobinemia/uria)
- + No agglutination & Coomb's -
- + Ddx
  - + Onion, propylene glycol (cats), water intoxication
  - + Heavy metal
  - + PFK/PK deficiency
  - + Hypophosphatemia
  - + Microangiopathy
  - + Infectious
  - + Massive envenomation, sprayed by skunk

## Immune-mediated hemolytic anemia (IMHA)

- + Extravascular (bilirubinemia/uria);
  - + Rarely intravascular or at level of bone marrow
- + Often autoagglutination & Coomb's +
- + Ddx
  - + "Triggered"
    - + Infectious disease
    - + Drugs
    - + Vaccinations
    - + Neoplasia
  - + Primary/autoimmune

Weis DJ, Tvedten H. Erythrocyte disorders. In: Willard MD, Tvedten H, ed. Small Animal Clinical Diagnosis by Laboratory Methods. St. Louis MO: Elsevier, 2012: 48-55.

# Missy

- 13-year-old fs Shih Tzu
- Occasionally collapses with exercise



# Triage

## + Quick PE



- + TPR = 101 F, 120 bpm, 20 bpm
- + BARH
- + Pink, moist MM, injected sclera
- + Normal pulses and chest auscultation

## + Assessment

- + Unremarkable PE, currently stable dog
- + History of occasional exercise-induced collapsing

## + Diagnostics

- + 1<sup>st</sup> tier: PCV/TP, BG, electrolytes? Full CBC/chem/UA? ECG?
  - + Consider exercise the patient and perform pre/post exercise diagnostics such as BG, lactate, K+, BP, ECG



Confirm Hct!  
PCV/TP = 92%/7.6 g/dL

# Focus

- + What's the client complaint? Collapsing
  - + What's the main problem? Erythrocytosis, likely polycythemia
  - + What's the emergency? Poor perfusion
  - + Do these all fit together; are we addressing all issues?
- 
- + What now?
    - + Restore perfusion > phlebotomy
    - + 4th tier diagnostics: thoracic radiographs/met check, abdominal imaging, +/- blood gas, erythropoietin concentration, echo



[Sludge Vectors by Vecteezy](https://www.vecteezy.com/free-vector/sludge)

# Treatment of erythrocytosis

## + Relative

- + Fluid therapy
- + Treat the cause of dehydration

## + Absolute

- + ER treatment = phlebotomy
  - + Dogs: 20 ml/kg, replace with same volume fluid, goal is to get PCV  $\leq$  65%, may need to repeat
  - + Cats: 10 ml/kg, same protocol, target PCV  $\leq$  about 50%
- + All cases may need long-term, occasional phlebotomy
- + IDEALLY, treat the underlying disease if there is a treatment
  - + Polycythemia vera: hydroxyurea
  - + Renal tumors > surgery if an option
  - + Other cancers, appropriate medical and/or surgical management as guided by an oncologist

Giger U: Polycythemia—Diagnostics and Management of a High Hematocrit.  
World Small Animal Veterinary Association Congress Proceedings 2018.

[Patient management >](#)

Canine | Shih Tzu | Female | 13 y

2019 **Dec 19** Dec 3 Nov 22 **Nov 21** 2011 Jul 19[Result Details v](#) **Hematology**12/19/19  
3:46 PM 12/3/19  
3:57 PM 11/22/19  
3:38 PM 11/21/19  
6:09 PM 

<b>RBC</b>	4.96	11.98	14.54	15.52	5.65 - 8.87 M/ $\mu$ L
<b>Hematocrit</b>	32.1	*71.8	*84.8	*91.0	37.3 - 61.7 %
<b>Hemoglobin</b>	9.0	22.3	*26.8	*28.5	13.1 - 20.5 g/dL
<b>MCV</b>	64.7	59.9	58.3	58.6	61.6 - 73.5 fL
<b>MCH</b>	18.1	18.6	18.4	18.4	21.2 - 25.9 pg
<b>MCHC</b>	28.0	31.1	31.6	31.3	32.0 - 37.9 g/dL
<b>RDW</b>	22.8	26.3	29.2	29.9	13.6 - 21.7 %
% Reticulocytes	0.1	0.1	0.8	0.8	%
<b>Reticulocytes</b>	6.0	16.8	119.2	118.0	10.0 - 110.0 K/ $\mu$ L
<b>Reticulocyte Hemoglobin</b>	22.9	20.7	18.5	18.6	22.3 - 29.6 pg
<b>WBC</b>	3.44	4.60	8.92	7.96	5.05 - 16.76 K/ $\mu$ L

# Henry

- 5-year-old m/c German Shorthaired Pointer
- Collapsing and heavy breathing with exercise
- Occasional cough



# Triage

## + Quick PE

- + TPR = 101 F, 180 bpm, 48bpm
- + Irregularly irregular heart rate, 2/6 systolic murmur
- + Pale to muddy, moist MM
- + RR and effort increase with minimal stress
- + Mild increase in bronchovesicular lung sounds
- + Moderate pulses

## + Quick assessment

- + Tachyarrhythmia
- + Probable lower airway issues, CHF?
- + Anemia?
- + Bordering not stable

# What diagnostics?

## + FIRST STABILIZE!!!

- + Oxygen
- + Consider light sedation
- + Consider diuretic such as furosemide

+ 1<sup>st</sup> tier diagnostics: ECG!!!, pulse oximetry, SBP, +/- PCV/TP/BG (CBC/chemistry/UA)

+ 2<sup>nd</sup> tier diagnostics: Thoracic rads or TFAST


+ 3<sup>rd</sup> tier diagnostics: Echo, proBNP, +/- refer to cardiologist

ECG = A-fib

Pulse ox = 93%

PCV/TP= 42%/6.8 g/dL


SBP = 90 mmHg



# Focus



- + What's the client complaint? collapsing, coughing
- + What's the main problem? Cardiac disease, +/- lower airway issue (including CHF)
- + What's the emergency? Poor perfusion, possibly oxygenation
- + Does this all fit together/are we addressing all issues/concerns?
  
- + What next?







# Echo



# Treatment

- + Oxygen
- + Diuretic: e.g., furosemide
- + Negative chronotrope: e.g., diltiazem, beta blocker, digoxin
- + Positive inotrope: e.g., pimobendan, digoxin
- + Longer term treatment: +/- decrease afterload, e.g., ACEI, ARB



## Follow-up?

- + Monitor:
  - + Hydration
  - + Renal chemistries and electrolytes
  - + SBP, ECG, RR & effort, rads, echo, +/- proBNP
- + \*\* Check for grain-free diet!!!\*\*

# Brutus

- 10-year-old MC Rottweiler
- Snoring more than usual for months, progressively worse
- Anxious this am and collapsed



# Triage

## + Quick PE

- + TPR 103.5 F, 130 bpm, marked stertor
- + Paradoxical chest movement
- + Muddy/greyish MM, moist
- + Normal pulses, marked referred upper airway noise
- + Panicked b/c can barely breathe

## + Assessment

- + Unstable, can barely breathe, panics with any handling
- + Cardiovascularly stable for the moment

# Diagnostics

## + STABILIZE THE PATIENT

- + Sedation
- + Rest, minimal stress
- + +/- Tracheostomy?
- + +/- Oxygen


- + 1<sup>st</sup> tier: Sedated oral exam – **\*\*BE PREPARED\*\***, pulse oximetry
- + 2<sup>nd</sup> tier: ?
- + 3<sup>rd</sup> tier: ?

# Sedated upper airway exam!!!



# Focus

- + What's the client complaint? Can't breathe
  - + What's the main problem? Soft palate mass obstructing airway
  - + What's the emergency? At risk of respiratory arrest
  - + Do these all fit together; are we addressing all issues?
- 
- + What now?
    - + SECURE AN AIRWAY!!!



# Stridor - Respiratory – upper airway



# Max

- 2-year-old MI Standard Poodle
- Pre-neuter labs
- No complaints



## Chemistry

5:18 AM

Glucose	61	72 - 175 mg/dL	
Creatinine	1.7	0.9 - 2.5 mg/dL	
BUN	27	16 - 37 mg/dL	
BUN: Creatinine Ratio	15.9		
Phosphorus	4.2	2.9 - 6.3 mg/dL	
Calcium	0.7	8.6 - 10.6 mg/dL	
Sodium	153	147 - 157 mmol/L	
Potassium	11.5	3.7 - 5.2 mmol/L	
Na: K Ratio	13	29 - 42	
Chloride	120	114 - 126 mmol/L	
TCO2 (Bicarbonate)	13	12 - 22 mmol/L	
Anion Gap	32	12 - 25 mmol/L	
Total Protein	6.4	6.3 - 8.8 g/dL	
Albumin	2.4	2.6 - 3.9 g/dL	
Globulin	4.0	3.0 - 5.9 g/dL	

# Triage

- + Quick PE

- + Completely normal

- + Assessment

- + Very abnormal, scary labs

- + Clinically stable patient ?

# Diagnostics

## + STABILIZE THE PATIENT

+ ?

## + Next diagnostics

+ 1<sup>st</sup> tier: Confirm results, NEW blood sample

+ 2<sup>nd</sup> tier: ?

+ 3<sup>rd</sup> tier: ?

# Focus

- + What's the client complaint? Routine surgery
- + What's the main problem? Potentially something major
- + What's the emergency? If real, complications of profound hypocalcemia and hyperkalemia
- + Do these all fit together; are we addressing all issues?
  
- + What now?

Chemistry		5:18 AM	
Glucose	61	72 - 175 mg/dL	
Creatinine	1.7	0.9 - 2.5 mg/dL	
BUN	27	16 - 37 mg/dL	
BUN: Creatinine Ratio	15.9		
Phosphorus	4.2	2.9 - 6.3 mg/dL	
Calcium	0.7	8.6 - 10.6 mg/dL	
Sodium	153	147 - 157 mmol/L	
Potassium	11.5	3.7 - 5.2 mmol/L	
Na: K Ratio	13	29 - 42	
Chloride	120	114 - 126 mmol/L	
TCO2 (Bicarbonate)	13	12 - 22 mmol/L	
Anion Gap	32	12 - 25 mmol/L	
Total Protein	6.4	6.3 - 8.8 g/dL	
Albumin	2.4	2.6 - 3.9 g/dL	
Globulin	4.0	3.0 - 5.9 g/dL	

EDTA  
Contamination

# Precious

- 8-year-old FS Westie
- Referred for sedated ear cleaning
- History of chronic allergies and skin issues
- “FYI, collapse and quick recovery recently
- Generally healthy, “normal lab work”
- You note hypoglycemia on RDVM pre-anesthetic labs

































2024 **Jan 30** Jan 23


Result Details [▼](#)



### Chemistry

1/23/24  
6:00 AM 

 	<b>Glucose</b>	<b>37</b>	63 - 114 mg/dL	
 	IDEXX SDMA	10	0 - 14 µg/dL	
 	Creatinine	0.7	0.5 - 1.5 mg/dL	
 	BUN	12	9 - 31 mg/dL	
	BUN: Creatinine Ratio	17.1		
 	Total Protein	6.7	5.5 - 7.5 g/dL	
 	Albumin	3.5	2.7 - 3.9 g/dL	
 	Globulin	3.2	2.4 - 4.0 g/dL	
 	Albumin: Globulin Ratio	1.1	0.7 - 1.5	
 	ALT	59	18 - 121 U/L	
 	ALP	125	5 - 160 U/L	



# Triage

## + Quick PE

- + TPR 101 F/120 bpm/20 bpm
- + BARH, downright rambunctious
- + Pink, moist MM
- + Normal pulses and chest auscultation

## + Assessment

- + Completely stable vs. on the edge of crashing



# Diagnostics

## + STABILIZE THE PATIENT

- + Does this patient require treatment for profound hypoglycemia?

## + Next diagnostics

- + 1<sup>st</sup> tier: Repeat/confirm BG on glucometer = 27 mg/dL
- + 2<sup>nd</sup> tier: Consider cortisol, bile acids, insulin panel, thoracic rads, abdominal imaging
- + 3<sup>rd</sup> tier: Thoracic rads, advanced abdominal imaging

**Switching gears...**



# Focus

- + What's the client complaint? Referral for sedated ear cleaning
- + What's the main problem? Possibly hypoglycemic
- + What's the emergency? None at the moment; risk of seizures and other complications of profound hypoglycemia
- + Do these all fit together; are we addressing all issues?
- + What now?
  - + Treat hypoglycemia?: diet, +/- drugs, +/- glucose CRI?
  - + Ddx for asymptomatic, profound hypoglycemia?



## Chemistry

1/23/24  
6:00 AM



	Glucose	37	63 - 114 mg/dL	
	IDEXX SDMA	10	0 - 14 µg/dL	
	Creatinine	0.7	0.5 - 1.5 mg/dL	
	BUN	12	9 - 31 mg/dL	
	BUN: Creatinine Ratio	17.1		
	Total Protein	6.7	5.5 - 7.5 g/dL	
	Albumin	3.5	2.7 - 3.9 g/dL	
	Globulin	3.2	2.4 - 4.0 g/dL	
	Albumin: Globulin Ratio	1.1	0.7 - 1.5	
	ALT	59	18 - 121 U/L	



## Endocrinology

1/23/24  
6:00 AM



Cortisol

2.8

2.0 - 6.0 µg/dL



Cortisol

2.8

2.0 - 6.0 µg/dL



2024 Jan 30 Jan 23

Result Details 



**Diagnostic result printed** 3/18/2024 10:25 PM

[Details >](#)



## Endocrinology

1/30/24

1:35 AM



**Glucose**

**36.0**

63.0 - 114.0 mg/dL



**Insulin**

**>200.0**


5.2 - 41.5 uIU/mL



**Insulin: Glucose  
Ratio**

**725**

14 - 43 RATIO



Pet Name:	<input type="text"/>	Appt. Date:	02/15/2024 01:00 PM
Species:	Canine	Report Date:	02/15/2024 - Update: 02/15/2024
Breed:	West Highland White Terrier	Presenting Complaint:	
Sex:	Female - Spayed	Met check, insulinoma	
Weight:	17.6 lbs/ 8 kg		
Patient DOB:	01/29/2014 (10 years old)		

#### Abdominal Ultrasound Report

#### Interpretation:

**Liver:** nodule may be benign nodular regeneration, but metastatic neoplasia is not fully ruled out.

**Kidneys:** consistent with chronic degenerative kidney disease with dystrophic mineralization vs. urolithiasis.

**Pancreas:** nodule most likely insulinoma given high insulin:glucose ratio, although other neoplasia (carcinoma) vs. nodular hyperplasia are not fully ruled out.

**Lymph nodes:** consider inflammation vs. reactive lymphoid hyperplasia vs. metastatic neoplasia.

**Bladder:** Normal wall, two discrete stones (8.3mm, 4.3 mm) as well as a few small calculi near the bladder neck.

#### Interpretation:


**Liver:** nodule may be benign nodular regeneration, but metastatic neoplasia is not fully ruled out.

**Kidneys:** consistent with chronic degenerative kidney disease with dystrophic mineralization vs. urolithiasis.

**Pancreas:** nodule most likely insulinoma given high insulin:glucose ratio, although other neoplasia (carcinoma) vs. nodular hyperplasia are not fully ruled out.

**Lymph nodes:** consider inflammation vs. reactive lymphoid hyperplasia vs. metastatic neoplasia.

**Bladder:** stones as noted, consider calcium oxalate vs. struvite vs. other.



## Clinical History:

Insulinoma, pancreas, all visible tumor excised

## Pathologist's Report

### INTERPRETATION:

Islet cell neoplasm, histologically grade:malignant

Mitotic count (per 2.37 sq mm): 4

Histologic tumor-free margins: Clear; the nearest peripheral margin is only 0.11 mm.

Vascular invasion: Not present

cell tumors that can hypersecrete insulin, and can cause hypoglycemia.

These tumors can be benign or malignant, depending on the infiltrative nature of the neoplasm. Most insulinomas in domestic animals are malignant. Circumscription and encapsulation are not reliable indicators of benign behavior in insulinoma, as circumscribed and solitary carcinomas have been noted to have metastasized into adjacent lymph nodes or commonly into the liver. Carcinomas tend to be larger, have a multilobular appearance, contain foci of hemorrhage and necrosis, have less uniform size and shape, and demonstrate vascular and lymphatic invasion (Maxie. (ed.). 2007. Pathology of Domestic Animals. 5th edition. Vol 2. pp. 421-423.).

This one is multicentric within the pancreas, and infiltrative.

Immunohistochemistry, islet cell tumors have been noted to stain with insulin (18/18) most consistently, with variable staining with other peptides produced by normal islet cells: somatostatin (14/18), glucagon (9/18), and gastrin (1/18) (Hawkins et al. 1987. Immunocytochemistry of normal pancreatic islets and spontaneous islet cell tumors in dogs. Vet Pathol. 24:170.)

# Take Home

- + Symptoms are often non-specific
  - + Your choice of diagnostics in the heat of an emergency should be based on science AND art
- + Do not skip the basics
- + If it doesn't make sense, slow down before you act
- + Think forward





**VMX**  
2025  
VETERINARY MEETING & EXPO

*We thank you,  
with all our **heARTS***

**Please rate your speaker and session in the app!**

Presented By:

**NAV**  
YOUR VETERINARY COMMUNITY