Bovine Virus Diarrhea Antigen Test Kit

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Bovine Virus Diarrhea Antigen Test Kit

For veterinary use only.

Rapid in vitro diagnostic for detection of bovine viral diarrhea virus in persistently-infected animals (BVDV PI-animals) using ear-notch (ear-tissue biopsy) and serum samples.

Precautions and Warnings

- Handle all biological materials as though capable of transmitting BVD virus.
- Do not pipette by mouth.
- Do not eat, drink or smoke where specimens or kit components are being handled.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Do not mix components from kits with different serial numbers.
- Do not use a SNAP® device that has been activated prior to the addition of sample.
- Do not reuse pipettes that have been in contact with a specimen.
- Do NOT pipette conjugate from its bottle using a pipette that has been in contact with a specimen.
- Betadine cannot be used as a disinfectant when obtaining ear notches, as this will interfere with the assay.
- Do not perform assay in direct sunlight.
- Do NOT cross-contaminate sample collection device(s). Soak in disinfectant before and after collection.
- ALWAYS use one hand to stabilize SNAP® BVDV device during activation step, to prevent splashing.

Storage

- Store the SNAP® BVDV (device and conjugate) at 2–8°C.
- All components, including specimens, must be at 18–30°C before running the test.
- Remove from refrigeration only the number of devices and conjugate to be used.

Kit Components

	Large Ear Notch (when ear notch is ≥ 1 cm)	Small Ear Notch (when ear notch is ≥ 2–3 mm)	Serum
Provided	1 or 5 bottles of anti-BVDV: HRPO conjugate (17 mL each).		
	10 or 50 SNAP® BVDV devices		
	Reagents contained in each device: 0.6 mL substrate solution • 0.4 mL wash solution		
Provided (separate accessory pack)	Capped ear-notch tubes	Sample tubes with attached caps	
	Transfer pipettes (for conjugate)	Pipettes (for conjugate)	Pipettes (for sample and conjugate)
Self provided	Not applicable		A heating device is needed (such as a heater block, water bath, etc.) capable of maintaining a temperature of 45±5°C.



NOTE: see table on page 33 for the description of international symbols used on the kit labels.

Sample Information

REMINDER: All samples must be at 18–30°C before beginning the test procedure.

	Large Ear Notch (when ear notch is ≥ 1 cm)	Small Ear Notch (when ear notch is ≥ 2–3 mm)	Serum Samples
Source	Ear notches may be obtained from animals of any age. Ear notch must measure at least 1 cm on one side.	Ear notch must measure at least 2–3 mm in diameter. (Typically used in Europe)	150 µL of serum is required. Serum should be obtained from precolostral newborn calves or calves older than 3 months of age to prevent the presence of maternal anti-BVDV antibodies from interfering with the SNAP* BVDV. [†]
Storage	Short Term: Ear notches can be stored dry at 2–8°C. If refrigerated storage is not available, ear notches can be stored dry at 18–30°C for 3 days. Long Term: If prolonged storage is necessary, the ear notch can be stored frozen at or below -20°C.	Short Term: Serum can be stored at 2–30°C for up to 7 days. Long Term: Serum can be stored at or below -20°C. All previously frozen serum specimens must be thoroughly mixed prior to removal of the volume to be assayed.	

Test Procedure (Large and Small Ear-Notch)


	Large Ear Notch (when ear notch is ≥ 1 cm)	Small Ear Notch (when ear notch is ≥ 2–3 mm)
1	If stored in a refrigerator, allow all components to equilibrate at 18–30°C for at least 30 minutes.	
2	Place the SNAP* BVDV device on a flat surface. The SNAP* BVDV device must be kept horizontal throughout the test to ensure accurate results.	
3	<p>Dispense: 1.5 mL of conjugate into the tube containing the ear-notch.</p> 	<p>Dispense: 450 µL of conjugate into the sample tube containing the ear notch. (To aspirate 450 µL, draw conjugate to the mark on the dropper provided in the Serum and/or Small Ear Notch Accessory Pack.)</p> 
4	<p>Swirl: Soak the ear notch in the conjugate.</p> 	

5**Incubate:**

At 18–30°C for at least **10 minutes** (range: 10 minutes–3 days).

**6****Pour:**


Pour or transfer the contents of the tube into the sample well. Some conjugate will remain in the tube containing the ear-notch. (Basin should be full, ~1 mL).

**Incubate:**

At 18–30°C for **1 hour** (range: 1 hour–3 days).

**Pour:**

Pour or transfer the contents of the sample tube into the sample well. Some conjugate may remain in the sample tube containing the small ear notch. If the small ear notch falls into the sample well, carefully remove notch prior to activating device.

**7****Watch:**

Sample will flow across the result window, reaching the activation circle in approximately **15–60 seconds**.


NOTE: Some sample may remain in the well.

Watch carefully for the sample or blue color to appear in the activation circle.

8**SNAP!**

When color FIRST appears in activation circle, push activator down firmly until it is flush with the device body. Push hard!

NOTE: Some sample may not flow to the activation circle within 60 seconds; therefore, the activation circle may not change color. In this case, press the activator if the sample crossed the result window. If the sample has not crossed the result window, repeat testing with another device.

9**Read:**

Read test result at **10 minutes** (range: 10–15 minutes).

NOTE: Color in the area of the positive control spot may develop sooner, but results are not complete until 10 minutes after the device has been activated.

Test Procedure (Serum)

Serum

1

If stored in a refrigerator, allow all components to equilibrate at 18–30°C for at least 30 minutes.

2

Place the SNAP* BVDV device on a flat surface. The SNAP* BVDV device must be kept horizontal throughout the test to ensure accurate results.

3

Dispense:

Dispense 150 µL (or 3 drops using a sample pipette held in a vertical position) of sample into a fresh sample tube.

4

Dispense:

Dispense 200 µL (or 6 drops using a fresh sample pipette held in a vertical position) of conjugate into the sample tube.

5

Mix:

Cap the sample tube and mix thoroughly by inverting it 3–5 times.

6

Incubate:


Place the capped tube of sample-conjugate mixture into a 45°C (range: 40–50°C) incubation device for **10 minutes** (range: 10 to 20 minutes).

7

Pour:

Pour the contents of the sample tube into the sample well.


8 **Watch:**

Sample will flow across the result window, reaching the activation circle in approximately **15–60 seconds**.


NOTE: Some sample may remain in the well.

Watch carefully for the sample or blue color to appear in the activation circle.

9 SNAP!

When color FIRST appears in activation circle, push activator down firmly until it is flush with the device body. Push hard!

NOTE: Some sample may not flow to the activation circle within 60 seconds; therefore, the activation circle may not change color. In this case, press the activator if the sample crossed the result window. If the sample has not crossed the result window, repeat testing with another device.

10 Read:


Read test result at **10 minutes** (range: 10–15 minutes).

NOTE: Color in the area of the positive control spot may develop sooner, but results are not complete until 10 minutes after the device has been activated.

Interpreting the Test Results

Spot Locations of Device when Removed from Bag


Assay Validity Requirements

- Positive Control Spot develops a blue color.
- Negative Control Spot shows no additional color development greater than the background.

Assay Interpretation

- Positive Results (PI Status)**—the sample is considered positive if the blue color that develops in the sample spot is **more intense** than the background.

Note: A positive result in the SNAP BVDV test, or other test methods including PCR, may not always indicate persistently infected animals. A definitive diagnosis requires a second test using a fresh sample taken three weeks later. This second sample should also be positive in order to confirm PI status.

- Negative Results (Non-PI Status)**—the sample is considered negative if the blue color in the sample spot is **equal to or less intense** than the background.

Invalid Results

- Background**—If the sample is allowed to flow past the Activation Circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No Color Development (Positive Control spot)**—If positive control does not develop color, repeat the test.
- Color Development (Negative Control spot)**—If the negative control develops color that is more intense than the surrounding matrix background in the device, repeat the test.

Attention Customers in the State of California, USA

Bovine Viral Diarrhea Virus (BVDV) is a reportable disease in the State of California. For all BVDV testing performed in California, you are required to submit a Monthly Summary Report to the California Department of Food and Agriculture, Animal Health Branch Headquarters, declaring any Positive BVDV Test results.

Please contact the **CDFA** at cavet@cdfa.ca.gov or see the List of Reportable Diseases at: http://www.cdfa.ca.gov/ahfss/Animal_Health/Disease.html for more details.

Sensitivity and Specificity

Relative Sensitivity and Specificity of SNAP® BVDV Using Characterized Positive and Negative Samples							
Study ID SIK-0384	Sample Type	Sample Size SNAP BVDV/Reference				Relative Sensitivity and Specificity 95% Confidence Limits	Kappa Statistic
		+/-	-/+	+/-	-/-		
		Total					
Small Ear Notch	106	0	0	130	236	Sen. 100% (95% CL 96.5%–100%) Spec. 100% (95% CL 97.1%–100%)	1.00
Large Ear Notch	209	2	1	204	416	Sen. 99.1% (95% CL 96.3%–99.9%) Spec. 99.5% (95% CL 97.0%–100%)	0.99
Serum	211	9	0	206	426	Sen. 95.9% (95% CL 92.3%–97.9%) Spec. 100% (95% CL 97.7%–100%)	0.96

Samples characterized by PCR or virus isolation. Results based on single sampling; PI-status not confirmed.
Data on file at IDEXX Laboratories, Inc. Westbrook, Maine, USA.

References

1. Palfi V, Houe H, Philipsen J. Studies on the decline of bovine virus diarrhoea virus (BVDV) maternal antibodies and detectability of BVDV in persistently infected calves. *Acta Vet Scand.* 1993;34(1): 105–107.

IDEXX Technical Services:

IDEDEX Europe

Tel: 00800 727 43399

Fax: 00800 433 99329

idexx.com

IDEDEX USA & CANADA

Tel: 1 800 548 9997 or 1 207 556 4890

Fax: 1 800 328 5461 or 1 207 556 4826

idexx.com

IDEDEX South Africa

Tel: +2711 803 3001/2/3

Fax: 00800 433 99329

idexx.com

IDEDEX Australia & New Zealand

Tel: 1 800 655 978

idexx.com

U.S. Vet License No. 313

Product code 5A05.22

¹Maternal anti-BVD antibodies, which can be passed onto newborn calves in the first 24 hours of colostrum intake, may interfere with the test, producing false-negative results from immunotolerant, persistently infected (PI) calves. Because the level of maternal antibody decreases as the calf ages, this interference can be prevented by specifying age requirements for animals tested by this method.¹

*SNAP and Test with Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. and its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2012 IDEXX Laboratories, Inc. All rights reserved.

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Dados arquivados na IDEXX Laboratories, Inc. Westbrook, Maine, EUA.
© 2012 IDEXX Laboratories, Inc. Todos os direitos reservados.
IDEXX Europe B.V.