# Colilert<sup>\*</sup> 250

06-27397-03

**IDEXX**

**For Technical Support, please call:**

North/South America: 1 207 556 4496/1 800 321 0207

Europe: 00800 4339 9111

UK: +44 (0) 1638 676800

China: +86-21-61279528

Japan: +03 5301 6800

Australia: 1300 443 399

**IDEXX**

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

[idexx.com/water](http://idexx.com/water)

# Colilert\* 250 Test Kit

## Introduction


Colilert\* 250 simultaneously detects total coliforms and *E. coli* in drinking water. It is based on IDEXX's proprietary Defined Substrate Technology\*. When total coliforms metabolize the Colilert 250 DST\* nutrient indicator, ONPG, the sample turns yellow. When *E. coli* metabolize the Colilert 250 DST\* nutrient indicator, MUG, the sample also fluoresces. Colilert 250 is able to simultaneously detect these bacteria at 1 cfu/250 mL within 24 hours, even with as many as 5 million heterotrophic bacteria per 250 mL present.

## Storage

Store at 2–30°C away from light.

## Presence/Absence (P/A) Procedure for Noncarbonated Water Samples

1. Add contents of one pack to a 250 mL sample in a sterile, transparent, nonfluorescing vessel.
2. Cap vessel and shake.
3. Incubate at 36°C ± 2°C for 24 hours.
4. Read results according to Result Interpretation table below.

## Result Interpretation

Appearance	Result
Less yellow than the comparator <sup>1</sup>	Negative for total coliforms and <i>E. coli</i>
Yellow equal to or greater than the comparator	Positive for total coliforms
Yellow and fluorescence equal to or greater than the comparator	Positive for <i>E. coli</i>

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and toward the sample.
- Colilert 250 results are to be read after 24 hours of incubation.
- However, if the results are ambiguous to the analyst based on the initial reading, incubate up to an additional four hours (but not to exceed 28 hours total) to allow the color and/or fluorescence to intensify.
- Positives for both total coliforms and *E. coli* observed before 24 hours and negatives observed after 28 hours are also valid.
- In addition, laboratories may incubate samples for additional time (up to 28 hours total) for their convenience.

## Procedural Notes

- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a negative control can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- If a water sample has some background color, compare inoculated Colilert 250 sample to a control blank of the same water sample.
- Colilert 250 is a primary water test. Colilert 250 performance characteristics do not apply to samples altered by any preenrichment or concentration.
- Colilert 250 has not been validated for use with flavored bottled water.
- In samples with excessive chlorine, a blue flash may be seen when adding Colilert 250. If this is seen, consider the sample to be invalid, and discard sample.
- Aseptic technique should always be followed when using Colilert 250. Dispose of in accordance with Good Laboratory Practices.

## Quality Control Procedures

One of the following quality control procedures is recommended for each lot of Colilert 250:

- 1.A. IDEXX-QC Coliform and *E. coli*<sup>2</sup>: *Escherichia coli*, *Klebsiella variicola*<sup>3</sup>, and *Pseudomonas aeruginosa*.  
B. Quanti-Cult<sup>3</sup>: *Escherichia coli*, *Klebsiella pneumoniae* and *Pseudomonas aeruginosa*.  
C. Fill three sterile vessels with 250 mL of sterile non-buffered oxidant-free water and inoculate with a sterile loop of ATCC<sup>4</sup> strains, *Escherichia coli* ATCC 25922/WDCM 00013 or ATCC 11775/WDCM 00090, *Klebsiella variicola*<sup>3</sup> ATCC 31488/WDCM 00206 and *Pseudomonas aeruginosa* ATCC 10145/WDCM 00024 or ATCC 27853
2. Follow the P/A Procedure above.
3. Results should match the Result Interpretation table above.

1. IDEXX P/A Comparator, catalog #WP104

2. IDEXX-QC Coliform and *E. coli*, IDEXX Catalog #UN3373-WOC-TCEC

3. Quanti-Cult cultures, catalog #WKIT-1001

4. American Type Culture Collection, 1-800-638-6597 atcc.org.

‡. *Klebsiella pneumoniae* (ATCC 31488/WDCM 00206) has been renamed to *Klebsiella variicola*.

\*Colilert and Defined Substrate Technology, DST are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Quanti-Cult is a trademark or registered trademark of Remel Inc.

Patent information: [idexx.com/patents](http://idexx.com/patents).

© 2017 IDEXX Laboratories, Inc. All rights reserved.

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© 2017 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.