# Colilert<sup>\*</sup>

06-12999-09

**IDEXX**

IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2017

---

**IDEXX**

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA  
[idexx.com/water](http://idexx.com/water)

**For Technical Support, please call:**

North/South America: 1 207 556 4496/1 800 321 0207

Europe: 00800 4339 9111

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: 03 5301 6800

Australia: 1300 443 399

# Colilert\* Test Kit

## Introduction

Colilert\* simultaneously detects total coliforms and *E. coli* in water. It is based on IDEXX's proprietary Defined Substrate Technology\*. When total coliforms metabolize Colilert's DST\* nutrient-indicator, ONPG, the sample turns yellow. When *E. coli* metabolize Colilert's DST\* nutrient-indicator, MUG, the sample also fluoresces. Colilert can simultaneously detect these bacteria at 1 cfu/100 mL within 24 hours even with as many as 2 million heterotrophic bacteria per 100 mL present.

## Storage

Store at 2–30°C away from light.

## Presence/Absence (P/A) Procedure

1. Add contents of one pack to a 100 mL sample in a sterile, transparent, nonfluorescing vessel.
2. Cap vessel and shake.
3. Incubate at 35±0.5°C for 24 hours.
4. Read results according to Result Interpretation table below.

## Quanti-Tray\* Enumeration Procedure

1. Add contents of one pack to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until dissolved.
3. Pour sample/reagent mixture into a Quanti-Tray\* or Quanti-Tray\*/2000 and seal in an IDEXX Quanti-Tray\* Sealer.
4. Place the sealed tray in a 35±0.5°C incubator for 24 hours.
5. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

## Result Interpretation

Appearance	Result
Less yellow than the comparator <sup>1</sup>	Negative for total coliforms and <i>E. coli</i>
Yellow equal to or greater than the comparator	Positive for total coliforms
Yellow and fluorescence equal to or greater than the comparator	Positive for <i>E. coli</i>

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Colilert results are to be read after 24 hours of incubation.
- However, if the results are ambiguous to the analyst based on the initial reading, incubate up to an additional four hours (but not to exceed 28 hours total) to allow the color and/or fluorescence to intensify.
- Positives for both total coliforms and *E. coli* observed before 24 hours and negatives observed after 28 hours are also valid.
- In addition, laboratories may incubate samples for additional time (up to 28 hours total) for their convenience.

## Procedural Notes

- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures. For example, samples run in other countries are incubated at 36±2°C for 24–28 hours.
- Colilert can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*<sup>2</sup> MPN tables should be used to find Most Probable Numbers (MPNs).
- If a water sample has some background color, compare inoculated Colilert sample to a control blank of the same water sample.
- If sample dilutions are made, multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- Colilert is a primary water test. Colilert performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- In samples with excessive chlorine, a blue flash may be seen when adding Colilert. If this is seen, consider sample invalid and discontinue testing.
- Aseptic technique should always be followed when using Colilert. Dispose of in accordance with Good Laboratory Practices.

## Quality Control Procedures

1. One of the following quality control procedures is recommended for each lot of Colilert:
  - A. IDEXX-QC Coliform and *E. coli*<sup>3</sup>: *Escherichia coli*, *Klebsiella variicola*<sup>4</sup>, and *Pseudomonas aeruginosa*
  - B. Quanti-Cult<sup>SM</sup>: *Escherichia coli*, *Klebsiella pneumoniae* and *Pseudomonas aeruginosa*.
  - C. Fill three sterile vessels with 100 mL of sterile nonbuffered oxidant-free water and inoculate with a sterile loop of ATCC<sup>5</sup> strains, *Escherichia coli* ATCC 25922/WDCM 00013 or ATCC 11775/WDCM 00090, *Klebsiella variicola*<sup>4</sup> ATCC 31488/WDCM 00206 and *Pseudomonas aeruginosa* ATCC 10145/WDCM 00024 or ATCC 27853.
2. Follow the P/A Procedure or Quanti-Tray Enumeration Procedure above.
3. Results should match the Result Interpretation table above.

**NOTE:** IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at [idexx.com/water](http://idexx.com/water).

1. IDEXX P/A Comparator, catalog #WP104; Quanti-Tray Comparator #WQTC, or Quanti-Tray/2000 Comparator #WQT2KC

2. Eaton, AD, Clesceri, LS, Greenberg, AE, Rice, EN. *Standard Methods for the Examination of Water and Wastewater*. American Public Health Association, 2005. Washington, DC.

3. IDEXX-QC Coliform and *E. coli*—IDEXX Catalog #UN3373-WQC-TCEC

4. Quanti-Cult cultures—IDEXX catalog # WKIT-1001

5. American Type Culture Collection 1-800-638-6597 atcc.org

†. *Klebsiella pneumoniae* (ATCC 31488/WDCM 00206) has been renamed to *Klebsiella variicola*.

\*Colilert, Defined Substrate Technology, DST and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. Quanti-Cult is a trademark or registered trademark of Remel Inc.

Patent information: [idexx.com/patents](http://idexx.com/patents).

© 2019 IDEXX Laboratories, Inc. All rights reserved.

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© 2019 IDEXX Laboratories, Inc. Wszelkie prawa zastrzeżone.