# Legiolert<sup>\*</sup>

06-0005740-04

**IDEXX**

**For Technical Support, please call:**

North/South America: +1 207 556 4496/+1 800 321 0207

Europe: +00800 4339 9111

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: +81 422 71 5921

Australia: +1300 443 399

Italy: +39 055 208450

**IDEXX**

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

[idexx.com/water](http://idexx.com/water)

# Legiolert® Test Kit

## Introduction

Legiolert is a test that detects *Legionella pneumophila* in water samples. This test is based on a bacterial enzyme detection technology that signals the presence of *Legionella pneumophila* through utilization of a substrate present in the Legiolert reagent. *Legionella pneumophila* cells grow rapidly and reproduce using the rich supply of amino acids, vitamins and other nutrients present in the Legiolert reagent. Actively growing strains of *Legionella pneumophila* use the added substrate to produce a brown color indicator. Legiolert detects *Legionella pneumophila* at 1 organism in 100 mL within 7 days.

## Storage

Store at 2–25°C away from direct light and humidity.

## Potable Water Sample Preparation

### 10 mL Protocol

1. Add 90 mL of sterile diluent (e.g., dechlorinated water, deionized water, phosphate buffer, Butterfield's buffer, or 0.1% peptone) to a sterile vessel.
2. Add 10 mL water sample.
3. Bring sample to room temperature.
4. Add contents of Legiolert blister pack to sample.<sup>A,B</sup>
5. Shake until contents are dissolved. Sample may remain cloudy.

### 100 mL Protocol<sup>1</sup>

1. Collect 100 mL of sample in a sterile vessel.
2. Bring sample to room temperature.
3. Using a hardness dip strip (supplied as part of the Legiolert Supplement kit<sup>2</sup>), follow the manufacturer's instructions to determine if your sample has low or high hardness. Consider 0-2 pads positive ( $\leq 7$  CaCO<sub>3</sub> (°d)) as low; 3-4 pads positive ( $> 14$  CaCO<sub>3</sub> (°d)) as high.
4. Add corresponding supplement (supplied as part of the supplement kit) as follows: Briefly agitate the supplement vessel immediately prior to use, then for low hardness samples add 0.33 mL and for high hardness samples add 1.0 mL.
5. Add contents of Legiolert blister pack to sample.<sup>A,B</sup>
6. Shake until contents are dissolved. Sample may remain cloudy.

## Nonpotable Water Sample Preparation

### 0.1 mL Protocol

1. Add 100 mL of sterile diluent (e.g., dechlorinated water, deionized water, phosphate buffer, Butterfield's buffer, or 0.1% peptone) to a sterile vessel.
2. Add contents of Legiolert blister pack to vessel.<sup>A,B</sup>
3. Shake until contents are dissolved. Solution may remain cloudy.
4. Add 0.2 mL of Legiolert Pretreatment into a sterile microtube.<sup>3</sup>
5. Add 0.2 mL of nonpotable water sample to the same microtube and mix well.
6. Incubate the microtube at room temperature for 60 seconds ( $\pm 5$  seconds).
7. Mix the microtube contents and **immediately** transfer 0.2 mL to the vessel containing Legiolert and mix well.

### 1.0 mL Protocol<sup>1</sup>

1. Add 100 mL of sterile diluent (e.g., dechlorinated water, deionized water, phosphate buffer, Butterfield's buffer, or 0.1% peptone) to a sterile vessel.
2. Add contents of Legiolert blister pack to vessel.<sup>A,B</sup>
3. Shake until contents are dissolved. Solution may remain cloudy.
4. Add 2.0 mL of Legiolert Pretreatment into a sterile tube.<sup>3</sup>
5. Add 2.0 mL of nonpotable water sample to the same tube and mix well.
6. Incubate the tube at room temperature for 60 seconds ( $\pm 5$  seconds).
7. Mix the tube contents and **immediately** transfer 2.0 mL to the vessel containing Legiolert and mix well.

## Quanti-Tray® Enumeration Procedure

### For Quality Control (QC) and Proficiency Test (PT) samples only:

1. Add contents of Legiolert blister pack to sample.<sup>A,B</sup>
2. Shake until contents are dissolved. Sample may remain cloudy.

### For all test samples:

3. Pour sample/reagent mixture into a Quanti-Tray/Legiolert tray.<sup>C</sup> Tap or flick the tray to remove air bubbles.<sup>5</sup>
4. Immediately seal tray in an IDEXX Quanti-Tray Sealer PLUS.
5. Incubate the sealed trays:
  - i. For all potable water or QC and PT samples, incubate at  $39 \pm 0.5^\circ\text{C}$  for 7 days.
  - ii. For all nonpotable water incubate at  $37 \pm 0.5^\circ\text{C}$  for 7 days.
6. Quanti-Tray/Legiolert **must** be incubated paper side down with the wells facing up in a humidified environment (approximately 85% humidity). Stack trays in alternating directions for added stability.<sup>F</sup> Do not stack more than 10 trays.
7. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the Legiolert MPN table provided with the trays to obtain a Most Probable Number. Multiply results by the dilution factor for the final MPN, if applicable.

## Result Interpretation

Appearance	Result
Any brown color, greater than the negative control (with or without turbidity)	Positive for <i>Legionella pneumophila</i>
Any turbidity, greater than the negative control (with or without any brown color change)	Positive for <i>Legionella pneumophila</i>
No brown color change and no turbidity, greater than the negative control	Negative for <i>Legionella pneumophila</i>

- Legiolert results are definitive at 7 days. In addition, positives observed before 7 days and negatives observed after 7 days are also valid.
- Turbidity alone indicates growth of *Legionella pneumophila* when the indicator may not develop as quickly or strongly.

## Procedural Notes

- For comparison, an incubated sterile diluent blank containing Legiolert reagent (negative control) can be used when interpreting results.
- Legiolert is a primary water test. Legiolert performance characteristics do not apply to samples altered by pre-enrichment or concentration.
- Aseptic technique should always be followed when using Legiolert. Dispose of materials in accordance with Good Laboratory Practices.

## Quality Control Procedures

**Note:** For QC and PT samples, please use the Quanti-Tray Enumeration section above.

1. One of the following quality control procedures is recommended for each lot of Legiolert:
  - A. IDEXX-QC *Legionella pneumophila*<sup>4</sup>: *Legionella pneumophila*, *Enterococcus faecalis*.
  - B. For each of the American Type Culture Collection (ATCC)<sup>5</sup> or World Data Centre for Microorganisms (WDCM)<sup>6</sup> bacterial strains:
    - i. Target organisms: *Legionella pneumophila* ATCC 33152/WDCM 00107 and/or ATCC 33156/WDCM 00180; streak the culture onto BCYE and incubate at 35°C for 48–72 hours.
    - ii. Non-target organism: *Enterococcus faecalis* ATCC 29212/WDCM 00087; streak onto blood agar plate and incubate at 35°C for 18–24 hours.
    - iii. For each bacterial strain inoculate a labeled sterile vessel containing 100 mL of sterile diluent with 10<sup>6</sup> to 10<sup>8</sup> cfu for target organisms or 10<sup>3</sup> to 10<sup>6</sup> cfu for non-target organisms.
2. Follow the Quanti-Tray Enumeration Procedure above.
3. Results should match the Result Interpretation table above.

**Note:** IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at [idexx.com/water](http://idexx.com/water).

<sup>1</sup> The potable water 100 mL sample protocol and the nonpotable water 1.0 mL sample protocol have not been validated in North America.  
<sup>2</sup> Legiolert Supplement kit, WLGT-SUP, 98-0005745-00  
<sup>3</sup> Legiolert Pretreatment, WLGT-PRE, 98-0007740-00  
<sup>4</sup> IDEXX-QC *Legionella pneumophila*, UN3373-WQC-LP 98-0009287-00  
<sup>5</sup> ATCC American Type Culture Collection 1-800-638-6597 • [atcc.org](http://atcc.org)  
<sup>6</sup> WDCM World Data Centre for Microorganisms • [wdcm.org](http://wdcm.org)

<sup>\*</sup> Legiolert and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. Patent information: [idexx.com/patents](http://idexx.com/patents).

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© 2017 IDEXX Laboratories, Inc. All rights reserved.
© 2017 IDEXX Laboratories, Inc. Todos os direitos reservados.