



Test **with Confidence**™

Intermixing reagents

IDEXX ELISA test kits are designed to work as a set of reagents. However, generic components, such as substrate, wash, and stop solutions can be mixed within kits having the same part number. Other reagents such as kit controls, conjugates and plates have been optimized to work together for that kit lot/batch. Should you have any questions on generic components or what can be used across kits, please contact your LPD Technical Services representative.

IDEXX Laboratories, Inc.

One IDEXX Drive, Westbrook, Maine 04092 USA, idexx.com
1-207-556-0300 1-800-548-6733 f. 1-207-556-4346