

# ImageVue peripherals and accessories

**Features**

**Benefits**

**Dimensions**

**What's in the box?**

**System computer—all ImageVue\* systems**



Dell\* OptiPlex\* XE4 SFF PC

- Processor: 12th Gen Intel\* Core\* i5-12400
- RAM: 16 GB
- Hard drive: 1 TB NVMe solid state drive
- Windows\* 11 Pro 64-bit operating system

- Ready to connect to in-practice networks or the web
- Improved startup time

- Height: 11.42"
- Width: 3.7"
- Depth: 11.5"
- Weight: 11.57 lb.

- Computer
- Standard wired mouse
- Standard wired keyboard
- Power cable

**Monitor—ImageVue digital radiography (DR) systems**



Barco\* MDRC-2222 Touchscreen Monitor

- DICOM calibrated
- Front consistency sensor
- Recommended modalities: All digital images except digital mammography
- Aspect ratio: 16:9
- Display area: 18.7" x 10.5"

- Resolution: 2 MP (1920 x 1080)
- Pixel pitch: 0.248 mm
- Contrast ratio: 1000:1
- Touchscreen capability
- Wide viewing angle (178°)

- Landscape dimensions with stand:
  - Height: 15.25"–19.5"
  - Width: 20.25"
  - Depth: 8"
  - Weight: 12.35 lb.

- Stand (riser and base)
- Power cable
- Video cables
- Quick install guide